

## Product datasheet for **SC208144**

### MAP7D1 (NM\_018067) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	MAP7D1 (NM_018067) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	MAP7D1
Synonyms:	PARCC1; RPRC1
ACCN:	NM_018067
Insert Size:	621 bp
Insert Sequence:	>SC208144 3'UTR clone of NM_018067 The sequence shown below is from the reference sequence of NM_018067. The complete sequence of this clone may contain minor differences, such as SNPs. <b>Blue</b> =Stop Codon <b>Red</b> =Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CAGTCCCCGCAGGTCACAGAAGTCCTTAAAGAGGGTTTGCCTTGGATCCGGGCACAGTTGTGAGGGCTC
CTCTGCATCACCTACCAGGATGTCTGGAGGAGAAAAAGACAGAACAAGATGGAAGTGGCCTGGGCCCC
TGGGGGTGGGTCTCTCTGTTGTTTTAACTGACCTTATAGACTGATGTCTCTTTGGCCGGAGCCAG
ATCTGCCCTCAGTGCATTTCGTGTGCTCGCACGGCAGACATCCCTTCTCCCCCATAACACATATACA
CTCACAGCTCTCTGGCCTCTTCCCTTGGGGAGGGGCCACCTGTAGTATTTGCCTTGATTTGGTGGGGT
ACAGTGGATGTGAATACTGTAATAGCTTGTGCTCAGACTCCTCTGCGTGGAGAGGGTGGGTGCAGGAG
GCAGACCCTCCCCCAAAGCCCCCTGGGGAGATCTTCTCTCTATTTAACTGTAAGTGGGGGATC
CCAGGTCTGGGGATGGGGACACCTTGGGCCACAGGATACTGGTTGCTTCAGGGGTACCCATGCCCCCT
GCCCTCGCCTGGAATCAGTGTTACTGCATCTGATTAATGTCTCCAGAAATAAGAATAATTCTGCCAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



[View online »](#)

RefSeq: NM\_018067.5

Locus ID: 55700

MW: 22.9