

Product datasheet for SC208197

PRP19 (PRPF19) (NM_014502) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PRP19 (PRPF19) (NM_014502) Human 3' UTR Clone
Symbol:	PRP19
Synonyms:	hPSO4; NMP200; PRP19; PSO4; SNEV; UBOX4
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_014502
Insert Size:	649 bp
Insert Sequence:	>SC208197 3'UTR clone of NM_014502 The sequence shown below is from the reference sequence of NM_014502. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GACAGAAGCCTCAAGTTCTACAGCCTGTAGGCCCTGGCCCTTCTGATGGAAGCTGGGCTCATCTCAGT
AGAGGGGTAGAATTAGGGTTTGGGGGGGGGGTGGGGGAATCTATGGGGGAGGGGCTCTGTGGGGT
GGGACATTCACATCATTTCACTCTGGTCTGAGTGGTGGCCTGAGAACCATGGTGGCATGGACCACCTC
ATCCATGCAACTCCAGGCCCATGGGAACGGATGTGGAAGGAAGAACTGTCACCCTCTTAAGGCCCAGG
GTCGGAGCCAGGCCTCTCCCTTCCTGTCGTTCAATGGACGTGGTGGTGGTGTTCACACCCATTTT
GTTGCAGTTCCTGTGAGACAGGAGAGGCTGAGCCAAGGGAAGTGTGAAGGGATGGGAGGAGGGCTTG
TGCAGGGTTTTGTAAGCAGTGATCTAGTTTCATTAAGGAAAGAAACAATAACCATTAACCACTCCCCG
TGTCTGTCTGCACCAGGAGCACCTGGGACTGGGAAGGTCAAGGGGAGGGAGCACACACTGGGACACTGG
CTTCCGGGAAGCCCATCTTCCTTTCTTTTACAGCTCTTACCCTTTTTTTTTTTTTTAATTGCACAGC
AGAAATAAAACAAATCTGCAGATGAAA
ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_014502.5
Summary:	PSO4 is the human homolog of yeast Pso4, a gene essential for cell survival and DNA repair (Beck et al., 2008 [PubMed 18263876]).[supplied by OMIM, Sep 2008]
Locus ID:	27339
MW:	23.3