

## Product datasheet for **SC208233**

### ZBTB22 (NM\_005453) Human 3' UTR Clone

#### Product data:

**Product Type:** 3' UTR Clones  
**Product Name:** ZBTB22 (NM\_005453) Human 3' UTR Clone  
**Symbol:** ZBTB22  
**Synonyms:** BING1; fru; fruitless; ZBTB22A; ZNF297; ZNF297A  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pMirTarget (PS100062)  
**ACCN:** NM\_005453  
**Insert Size:** 624 bp  
**Insert Sequence:** >SC208233 3'UTR clone of NM\_005453  
The sequence shown below is from the reference sequence of NM\_005453. The complete sequence of this clone may contain minor differences, such as SNPs.  
**Blue**=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GGCTTCGGTGGAGGTGGAGGAGCAAACTGAAGGGGCAGGCTACTGGGGTGGGGTAGCTTTCGGGAAAGG
GAATAAGGAGCACGATGCAAGGGCGCTGTGGCCCCGGGTGATCTCCACCACACTTACTGTCTTCCTT
TATCTCTGTGGACTTGTATATATTCTGGAAGGGGAACCACAGTTTCACCATCGCCCGCCATTCTACTA
CTCAACCCCTCCCCCAAGGATTTCCAGAATAAACCCCTTCCCTTCCCTCTGATGGGTACACTGAAG
CCCTGCTCCACAGAGTAGATTGCACATGGAGGGAGGAGGGGCGTGTGAACATCCTGCAGTCAC
AGGGTCAGGGGTGAGGTGGTGTAGTCTGTGCCTGAAGTCTGTGTTGTGTTGTCGTTGGAGACAAGGCC
TTTGAGCCCCACCTTGTCTAGAACCTACCCCTCTCAAGGATGCGCTCTTTATTTCTACCCTGTCTC
TCCCCGCCACCCCGACTTCCCGTGAAATTTCCAACCTCGGTTCTCATGGAGGAGTGGGTGGAGACAAG
GAGGGAGTAAGTCGTAGGAGTACAAGGTTTTTATTTTTTTAACAGTGATAAAAATTTATTGGTCAT
TTA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

**Restriction Sites:** Sgfl-MluI  
**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

|                    |  |
|--------------------|--|
| <b>Components:</b> | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials. |
| <b>RefSeq:</b>     | <a href="#">NM_005453.5</a>  |
| <b>Summary:</b>    | May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]  |
| <b>Locus ID:</b>   | 9278   |
| <b>MW:</b>         | 22.7   |