

Product datasheet for SC208246

KIF9 (NM_001134878) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	KIF9 (NM_001134878) Human 3' UTR Clone
Symbol:	KIF9
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001134878
Insert Size:	684 bp
Insert Sequence:	>SC208246 3'UTR clone of NM_001134878

The sequence shown below is from the reference sequence of NM_001134878. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
ATGGGCCTCCAGCAGGCACATAGAAAATAGAACCTCATCGCCAGTACCTTAAAGGACAAGACCAGCAAC
TCCCACCTACTGTAGTGGAGCTGCTCAACCACCTGCCAGAGCTGCAGCCCCCTCTACTCCAATGCTGG
GACCCAGCACAGAGAACACATTTGGCCTGCATGTTGGGAGGAGCATCCTCCAAGGACAACCTTGCTCAT
CTCCACAGAGCACTTTGGGTTTTAATTCAGTGTCTTATATGCAGGGACAGGATAAAATAACTTTCTAGT
TTGGACTTTTCAGTGCTCATCTGAGAACTTTGTCTTATTTCTGTCTGAAAGGACATATTCTCTATGTAT
GTGTGAGGTCAGAACTGAGTCTCAAGTAGGAAGGCAGAAAAGTAGCTCAGCCAGGCTCCTCCTGGCAGG
CCTTTGTCACTGGTCAGCCCCCTTCTGGATAGACATCCTAGGGTGGGGTGGTGAGCTCGGTCTGA
GGAGCTCTGCTTCTACAGCTCCCTCAGACCCTAAGTTCTCTCAAGACATCTTTCTGTGTTTCTGAGC
CAGTGAAGAAGACAGGAGGAGGCCCTCCTTCATTTGGCAACCAGGGGAGGATAAGGTCTAGCCACAGAG
ACCAGAAAGGGATGAGGCAGAGCTGGTCCATGTACCAGGAATAAATGAGCAGCTTCTGGGA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001134878.3</u>
Locus ID:	64147
MW:	25.3