

Product datasheet for SC208248

C9orf85 (NM_182505) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C9orf85 (NM_182505) Human 3' UTR Clone
Symbol:	C9orf85
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_182505
Insert Size:	645 bp
Insert Sequence:	>SC208248 3'UTR clone of NM_182505 The sequence shown below is from the reference sequence of NM_182505. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GACACAGGAGGAGACCATCAAATGAATTAAATCACTGTATTAAGTCTGCCGGGCACAGTGGCTCAC
GCCTGTAATCCCAACACTTTGGGAGGCCAAGGAGGGTGGATCACCTGAGGTCAGGAGTTCGAGACCAGC
CTGGCCAACATGGCGGAACCCCATCTCCACTAAAAGTACAAAAAATTAGCTGGGCGTGGTGGCTCATGC
CTGTAATCCCAGCTACTCAGGAGGCTGAGGCAGGAGGATTGCTTGAACCTGGAGCGGAGATTGAAGT
GAGCTGAGTTTCGTGCCATTACACTCCAGCCTGGGTGACAGAGTGAGACTCTGTCTCAAAAAAATAAAA
TAAAAAGTCAATTTAGAATGTGAAATTCTGACCACCTTTTGGCTTTGAGTATTTCCAAAAGATATTTG
AAATCCTAATGAGGAAATCAGAAAAAGCTATGAAAAATAGACAAATTCATACATGAACAATATAAAT
TGTGTATATTACTTAACATCAAATAAACAAGATTCAGAATTGATGGTTGTATAAGAACTAGCTCATGT
AAAAATAAAATAACATTATTACATTGCCTCAAAAATTGGTCCTCGGTAAGTGCCTTTGATAAATGATC
TCAAATAAATGACTCATAGCAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_182505.5](#)

Locus ID: 138241

MW: 24.3