

Product datasheet for SC208290

KIAA1841 (NM_032506) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: KIAA1841 (NM_032506) Human 3' UTR Clone
Symbol: KIAA1841
Mammalian Cell Selection: Neomycin
Vector: pMirTarget (PS100062)
ACCN: NM_032506
Insert Size: 643 bp
Insert Sequence: >SC208290 3'UTR clone of NM_032506
 The sequence shown below is from the reference sequence of NM_032506. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GTGCTCTTGATGACAGTGTGCTGTGTAAACGAATTAATACATGAAAGGACTTGAATACTGCTTGAC
ATATCGTGAGGCAGAGACACAAAGATATGGAAGAACATTGGACAGATTACATCGATATATTAATTAAG
GTGATTTTCAGAGTTGGAGCAGGAGAGGAAGCCTTATAGATATTTTTTTCTTTTGCTTGCTTGATT
TTCTTTTAAATTATTATTCAGTATTTCAAGTAGACAGAAAAGTACAGAGGATAATACAACAAATATCCA
GGTACCCACCCAGCATTTTTACATAGTAGCATTTTGATGTTTTGAAATAAGACATAAACATTCCAA
ATACCCTATCACCGAGGATCCTGCATCCTTCTTTTCCAGAAAGTGAACATTTCTCTGAATTTGATGC
TTATTATGCCCATGCATGTTGTTATACTTTTACTACATATTTTATGTAGCTACAGCATAATATGCTATG
TTTACTTTATATAGTGACCTCATACTATACATATCCTCATGCAATTTTTTAGACTCTGTTTGAGAATTA
TATATTTTGATTAGGTAGCTTAAGTTTATTTAGCTACTCAGTAGTAATGTATTCTATTAATCTACTA
ATATATTTTCTGTCTATGGA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI
OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_032506.4</u>
Locus ID:	84542
MW:	25.5