

Product datasheet for **SC208308**

C11orf49 (NM_024113) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C11orf49 (NM_024113) Human 3' UTR Clone
Symbol:	C11orf49
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_024113
Insert Size:	641 bp
Insert Sequence:	>SC208308 3' UTR clone of NM_024113 The sequence shown below is from the reference sequence of NM_024113. The complete sequence of this clone may contain minor differences, such as SNPs. Red=Cloning site Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCAA**GCGATCGC**

CTCCACTGAAGAGACAGACGAGTCGGAGACT**TGA**GGAGTCCAAAGGGTCCTGCCACAGCGCCCTGTACC
 TGCTCCCACCCAGCCCTTGGTGTGCCACCCAGCCTCCTCTCCAGCACCTTGCTGTGCTGCCCTCTGCTG
 CTGACAAGGTGAATAACAGCCCCAAGACCAGCCAGAGGGGCTCTGATGATCAGCCAGCCAGTGGCCCCG
 GAAGGTGAATGGCTGCTCTCCCTGGCCCTATCAGCCTGTGAACCTTCACTTAGGCCCAAGCTGACAGAC
 TGTGCTGAGGCCACCTTGTACGCCGTAGCCTGTTAGTCCTCCTAACCTCTTAAGAGCAGTCTCTCTGA
 GCCAGCCTCTGCGGGTCCCCAATAAGGTTTCATCTCTCACAGCAACTCCATTAAGGGGGAGAACCCGAA
 TAGCCACGCAGGGCCTTGACCATCAAGGGTGACACCTGCGACGCAAGTACCAGGAGGACATAACCGCTG
 TGGCCTGTTGGAGAACAGCCAGTAGCCTTGTAATATGAAGGGTGGGCCAGAAGATGATTTCACTTGCAA
 AAAGTCTTCAAGTCTTGACCCCTTTGTGTCTAATAGCTAAACAAACATGTGAAACGAATAAAAAGTCCC
 TCATGTCTGGT

ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCG

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_024113.3</u>
Locus ID:	79096