

Product datasheet for **SC208334**

PTPN13 (NM_006264) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: PTPN13 (NM_006264) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: PTPN13
Synonyms: FAP-1; hPTP1E; PNP1; PTP-BAS; PTP-BL; PTP1E; PTPL1; PTPLE
ACCN: NM_006264
Insert Size: 640 bp
Insert Sequence: >SC208334 3'UTR clone of NM_006264
The sequence shown below is from the reference sequence of NM_006264. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CAAAAACAGCAGCCTCAGCTTCTGAAGTGCATGAAAAGAGCCTCTGGATGCATTTCCATTTCTCTCCT
TAACCTCCAGCAGACTCCTGCTCTCTATCCAAAATAAAGATCACAGAGCAGCAAGTTCATACAACATGC
ATGTTCTCCTCTATCTTAGAGGGGTATTCTTCTGAAAATAAAAAATATTGAAATGCTGTATTTTTACA
GCTACTTTAACCTATGATAATTATTTACAAAATTTAACACTAACCAAACAATGCAGATCTTAGGGATG
ATTAAAGGCAGCATTGATGATAGCAGACATTGTTACAAGGACATGGTGAAGTCTATTTTTAATGCACCA
ATCTTGTTTATAGCAAAAATGTTTTCCAATATTTAATAAAGTAGTTATTTTATAGGGGATACTTGAAA
CCAGTATTTAAGCTTTAATGACAGTAATATTGGCATAGAAAAAGTAGCAAATGTTTACTGTATCAAT
TTCTAATGTTTACTATATAGAATTTCTGTAATATATTTATATACTTTTTCATGAAAATGGAGTTATCA
GTTATCTGTTTACTGATCATCTGTTGTAATCATTATCTCACTTTGTAATAAAAAACACACCTTA
AAACATGAACAAGCCAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



[View online >](#)

RefSeq: [NM_006264.3](#)

Summary: The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This PTP is a large intracellular protein. It has a catalytic PTP domain at its C-terminus and two major structural domains: a region with five PDZ domains and a FERM domain that binds to plasma membrane and cytoskeletal elements. This PTP was found to interact with, and dephosphorylate, Fas receptor and I κ B α through the PDZ domains. This suggests it has a role in Fas mediated programmed cell death. This PTP was also shown to interact with GTPase-activating protein, and thus may function as a regulator of Rho signaling pathways. Four alternatively spliced transcript variants, which encode distinct proteins, have been reported. [provided by RefSeq, Oct 2008]

Locus ID: 5783

MW: 24.8