

Product datasheet for **SC208343**

UBQLNL (NM_145053) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	UBQLNL (NM_145053) Human 3' UTR Clone
Symbol:	UBQLNL
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_145053
Insert Size:	651 bp
Insert Sequence:	>SC208343 3'UTR clone of NM_145053 The sequence shown below is from the reference sequence of NM_145053. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
ACACAGATTTCTGATCTGCTTAGTGCTTAGGCAACCCTAAAGCATTCCAAGCAATATTGCAGATTGAGC
AGGCCCTCCAGCTGCTGGCCACAGAGGCTCCTGTTCTTCTGCCTTGGGTTGCACCCTACCTATGGGGCC
TGGGTTGGCTTCTGCCCCAGCTGCAGCTATCCTGACACAGTGCCCTGTTCTTCTGGAATGTTTCAGATA
CAGCTGAGCCCAAGGGACCTGAGTGCTGCCACAAGCCTGGAACAGTCCTGCAGAGGCTACAATCCCCGG
ATGGGGACCCCTTCCACCCTCTGCAAGCTCCTGAGATTTGTTTTAGCAAACAGATGGATTCTCTCCAGG
CCATGGGATTTGGAACCACTGCCAATCTACAGGCACTCATTGCTACTGAAGGGACACCAATGCTG
CTATCCGCAAGCTCAAGAGATCCCAGAGATTCTAACCACCATGCCTACTTGTGCTTGCTACCTGCCT
GCTGACCCACCTGACCATCTCATTTGCCTTTTGACCTTTCTGATGCTTCCAGCCAGGAGAAGTCCTG
GAATAAGAGTTATCAACCAATGTGTCTTGTACTGAATAATAGATCATTGGTCGTGGCTGAAACATCTGT
CAATAAAATGGCTACACTCACTTGCTGCGA
ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_145053.5</u>
Locus ID:	143630
MW:	24.3