

Product datasheet for **SC208344**

GLB1L3 (NM_001080407) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	GLB1L3 (NM_001080407) Human 3' UTR Clone
Symbol:	GLB1L3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001080407
Insert Size:	625 bp
Insert Sequence:	>SC208344 3' UTR clone of NM_001080407 The sequence shown below is from the reference sequence of NM_001080407. The complete sequence of this clone may contain minor differences, such as SNPs. Red=Cloning site Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCAA**CGCATCGC**

ATCTACAGACAAGCCACGCT**GTAA**AACTGTGTCTGAACATTTTTTTTTTTTTTTTGGAGATGGAGTCTCAC
 TTTGTCGCCCAGGCTGGAGTGCAGTGGCACAATCTCTGCTCACTGCAAGCTCAGCCTCTCGGGTTCACGC
 CATTCTCCTGCCTCAGCCTCCCCAGCAGCTGGGACTACAGGTGCACGCCACCACGCCTGGCTAATTTTTT
 GTATTTTGTAGTAGAGATGGGGTTTCACCAAGTTAGCCAGGATGGTCCCAATCTCCTGACCTTGATCTG
 CTCTCCTCAGCCTCCCAAAGTACTGGGATTACAGGCGTGAGCCACCACTCCCGGCCGTGAACATATTTTT
 TGGGTTGCTGGAGTTCATCTATAAGTCATTTTTGAGGAATAAGATTTATGTTAAGACTATCAAACACAGT
 GTTGCTACAATAGCAAAAATGTGAAAATAACAACAACAAAAACAGCAGAGGAATTGTTATGTATTTT
 GTAGTCTATCTATGATGCCTATTTTAGGCTTTAAAAAGTCTTCAAATCTTAATGACTGATTTATC
 TAGTTAAATGCTTAATCCTTAGCAGGCTCTTATCTTTAATTAACGTGCCTTTGAGTAGATGTG

ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCG

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001080407.2</u>
Locus ID:	112937