

Product datasheet for **SC208348**

TMEM125 (NM_144626) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: TMEM125 (NM_144626) Human 3' UTR Clone
Symbol: TMEM125
Mammalian Cell Selection: Neomycin
Vector: pMirTarget (PS100062)
ACCN: NM_144626
Insert Size: 650 bp
Insert Sequence: >SC208348 3'UTR clone of NM_144626
 The sequence shown below is from the reference sequence of NM_144626. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
ACCATATCCAGCATTGCCAGCCTCATCTGACGGAGCCAGAGCCGTCCTTCTCTCACAGCGGCCTCAGC
GTCCCCAGAGCCGAGCCAGGGTGTGAGTGCATGTGAACGTTGAGTACACATGAGTGCCTGTATGCCCCC
AGGCTGGGTGAGCTCTTCTGTGGATTGCATGGCGTGTGATTAAGGCCCATGTGTTCCACACATCCAC
ATCATGGGAAGGTTAATGTGTGCCTCCTTGGAACTGGGTGTTGGTGCCATGGAATTCCTCTCTGTAT
CTCAGGTCAGTAGGCGCAGAAACGCCTCATGATGAAGATTCTTGAGCCCCATTTCCAAGACCCCTCACA
TCCAATCCTGTCTGTAACATCCATCAAGATTTCATAGGGGTGACTGGTGCCACCCAAGACTGCAC
CAGTGCCTGCTCATTGAGGAGAGTAAGTCTGGCCAGGCAGAAAGAATATGGGCTCTGCAATGAGACAG
ACCTGGAGGGGACTCTCCGTTGAGCACTAGCAGCTGGAGGAGTTGGGAGTTCATGGCTATCATGGTTG
TGTTAATCGATTGTGGGGATGAAATGTCATTGTGTATGGAAGGCGGGGCTCATGGCTGATTGGCAATAA
AATGGCGGCTGCCGTTGTCATTGTCTCCA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_144626.3</u>
Locus ID:	128218
MW:	23.4