

Product datasheet for SC208384

MAP11 (NM_018275) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	MAP11 (NM_018275) Human 3' UTR Clone
Symbol:	MAP11
Synonyms:	C7orf43; MCPH25; TRAPPC14
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_018275
Insert Size:	618 bp
Insert Sequence:	>SC208384 3'UTR clone of NM_018275 The sequence shown below is from the reference sequence of NM_018275. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
AAGGTCCTGGTGGTGGAAACCGTCAAGTAGCACCCTGCCAGCTCTGTTCCCTCTTACACTCCAGAGACC
CAACGCCCCCAGAGGGGATCCCTGCTCCCGGGCTGTGCCTCCCTGGGATGCCTCCAGACGGGGGTGA
AGAGGCCTGGCAGAGCTGCCTGTCTGTGTCTGCTGATGAGGGATGGGGGAAGAAGCTGTGAAGTGGGC
GGGCATGGCTGGGACTAAGCCACCAAGTATCCCCGAGTTCCTGTGGGGGGGGCTGGCCCCACCCCTAGG
CCAGGGCAAGGGTCCCAGAGCTCCCTGTCCCCGGCCCTTTACCCTGGTTCTGAGTTTACAAAGTCTC
TTCTCTATTCCCCTTGAGTTCTTTCCCACCTCTGACATTCCCTCCCTCCCTCCCGCAGGGCTGAGATTA
GAGGGTGGTGTAGGCTAAGGGCCCTGACAGTGACCTTCTGTCTCAGGGGTTGGGGACAGGGCCAGGT
AGCCTCTGCCCTTATGTTTACGTTTGACGCTGAAGCACTTTAATTTTTTTTTTTTGGTCTGTTCC
TGTAACCTATTTTCCAATATTGCTTCCAAGTAAATAAGACTATTAAATGCCTGTTTCAGGAGGGA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_018275.5
Summary:	Microtubule associated protein wich seems to play a role in spindle dynamics and cell division (PubMed:30715179, PubMed:30447097). Modulates YAP1 activity as transcriptional regulator (PubMed:30447097).[UniProtKB/Swiss-Prot Function]
Locus ID:	55262
MW:	22