

Product datasheet for SC208392

KRT71 (NM_033448) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	KRT71 (NM_033448) Human 3' UTR Clone
Symbol:	KRT71
Synonyms:	HYPT13; K6IRS1; KRT6IRS; KRT6IRS1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_033448
Insert Size:	653 bp
Insert Sequence:	>SC208392 3'UTR clone of NM_033448 The sequence shown below is from the reference sequence of NM_033448. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
AGTGCACCCTCCAAGAAAACAGTCGGTAGAGAAGACTGCCCCGGGCCCCGCCTCATTCCATGACCCGG
CTCTGGATCCCACACTGTACTTCCCACAGCCCACTCTCAGCTCCATCTCCACCCTGCTGGTCTGCTCC
CATACACCTGGCACTGGCCTTGGCCACCACTTCTCCAGCCTGTGTCTTCTGATCCTGGGAAGGCCT
GGATGACCAAGCTTGGTAAATTCTCCCTGTACACACCCTATTAACCTTGGCTGTGGTCCCCCAGC
TACACCACCAGCCAGGTCCTGGCTGCCAGCTTTCCTCCTCTGCCCGGCCTCTAGCGCAGTCGCTAACT
ACTCTGTCTGGGCTCCCTGGGTCTCTGCCAAGGCCCGCACACTGGGACCTAGCATAGTTCCTGCCT
ATGCCAGGAGCTGGCTCTGTGTTTAAGAAAAGGAGGACTGAAGGACAAACAACCAAGAGTGGCCAGTC
CCCACCCCCACATCTAGCTCAGTCTCAATCTGAGTGGGACCAAGTGCAATTCAGGGCCTTTTCTCCA
CTCACCTGCACCCAGAAGCAGAGAAAAGCAGGCACTGTTCACTTTTCTTATTCTTAATGGCCTTCCT
CTGTTGCAACCTCAATAAACAGCACAATCTCA
ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_033448.3</u>
Summary:	Keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into epithelial keratins and hair keratins. This gene encodes a protein that is expressed in the inner root sheath of hair follicles. The type II keratins are clustered in a region of chromosome 12q13.[provided by RefSeq, Jun 2009]
Locus ID:	112802
MW:	23.2