

Product datasheet for **SC208393**

ZNF841 (NM_001136499) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF841 (NM_001136499) Human 3' UTR Clone
Symbol:	ZNF841
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001136499
Insert Size:	672 bp
Insert Sequence:	>SC208393 3'UTR clone of NM_001136499 The sequence shown below is from the reference sequence of NM_001136499. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GTTGTCCTAACCTCTGGGATCCCCAAATTATTTATACTTACTCATATAGCTTGTATTTGTCCTTTCC
CTTTGAAGTCTCATGTGGAATTGTAGTAATCTCCAGTATTGGAGGTGGAGCCTAGTGGGAGGTGATTGG
ATGATGGGGTGGCTTTCTCAAGAAAGTTTGTAGTACCATCCCCTTGTTGCTGCTTGTGATAGTTCTC
ATGAGGTCTGGCTCATGCCTGTAATCCCAGCACTTTGGGAGGCTGAGGTGGCGGATCACTTGAGGTCA
GGAGTTCAAGAGCAGCCTGGCCAACATGGTGAAACCCTGTTTCTGCCAAAAATACAAAAATTAGCTGGG
CTTCGTGGCGGGTACTGGTAGTACCAGCTACTCGAGAGGCTGAGTCAGGAGAATTGCTTGAACCCTGGA
GGCGGAGGTTGCAGTGAGCCGAGATCGCGCCACTGCACTCTAGCCTGGGTGACAGAGTGAGACTCCATC
ACGCACACGCACGCACGAAAGTGAATGGCACCTCCCTCTGTTTGTCTTCTGGCCATGTGACATGCCC
TTCACCTCCCACCATCATTTTAACCTTCCTCAGGTTTCCCGAGAGGGTGAACAGACGTCAGCGCCATGT
TTCCTATAAAGCCTGCAAAATTGTGAGTCAATTAACCTCTTTATACATTA
ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001136499.2
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	284371
MW:	25.5