

Product datasheet for SC208445

WDR20 (NM_181302) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	WDR20 (NM_181302) Human 3' UTR Clone
Symbol:	WDR20
Synonyms:	DMR; FLJ33659; MGC33177; MGC33183
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_181302
Insert Size:	676 bp
Insert Sequence:	>SC208445 3'UTR clone of NM_181302 The sequence shown below is from the reference sequence of NM_181302. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
CCTGGTAAAGTGGTAAGTTTAAATCCTAATGCTGCACCAGATCTAGAACTTGAATAGGTAGTGACTTT
TTTCTTTTTCGTGGGAGGGGTGGGTGTACAATGAATGTGAATGACACTTCTTATTCTTAATGTAAATC
TCAATGCATCAGAGCCATAATTTTGGATACTGCATGCCATGTAATTCTGAATCATTGATAATTTACCT
TAGAGCATTTAAAAAATAATCAAATAATTGCCAGCCAAGTCAGTCATCCTCCTGGGAGTATATAG
AGTCCCAAGTTAGCGCTCCTGTATTAGACTATTTCAATTTTAGGAAAATCATGACCATGTGGGGAAC
AATGACTTTAAAAATGCTGAAATTAATTTATGCTTTAACTGGAATATTTTTTGCTTAACTACTCAATT
AGAATATTGTACACCTGATCAATGTGTGTTTCAGCACAGATGGCCATGAATTGTCATTTATAGTCCAATT
TTTTATCTTAATCATAAAATGTTTAGGAATCTATGAAATTTAACTTTAGGAACAAAACGTTTAGCAGGG
TTGATTGATATTATTTTACATTGTTCTGGCAATCCACAGAAAGAGAAGAGCCTTAATTTTTAAACCC
ATTTTAGTCATTTTATGACAATTAAGTTGTTTAATAAACATCTTTTTTCAAAGA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_181302.1</u>
Summary:	This gene encodes a WD repeat-containing protein that functions to preserve and regulate the activity of the USP12-UAF1 deubiquitinating enzyme complex. Multiple alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jun 2011]
Locus ID:	91833
MW:	26.4