

Product datasheet for SC208449

ARHGEF38 (NM_017700) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ARHGEF38 (NM_017700) Human 3' UTR Clone
Symbol:	ARHGEF38
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_017700
Insert Size:	661 bp
Insert Sequence:	>SC208449 3'UTR clone of NM_017700 The sequence shown below is from the reference sequence of NM_017700. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TTGAGCCACTGTATCCAGTCCTTAAAGTAAAGGCCCTTTCAAATGATGATTCCCATCTCCTCTCAGTTGC
CTAGCAGGGAACATTTTAAATGGATGTAGATGAAAGGTCTCACATAAATCCTATGTTTTATGAGACTTG
CTGGGAGCTCTGCTTTGCATTCCCTTTATAAAAGCTGACATGCCAGAAGCCCTGATTGACTTTTTTTC
CCCCTGCGAGAATGACTAAAAATAACATGGAAGAAGATTTAGAGCTCTGCAGCGATTGAAAAATGCAAT
ATCAAAATATAAAATGTGGAAGAAAAGCCTCTTCTTAAAGCTATTGTAACCTTGCCTGGCCCCACGTAGT
TCAAGGATTATGTGAGATAACACGTGGCCCCATGACCACTGGAGCACATGGGTTAATGGAGTTAGGGGA
ATGGCCTACAACCTCTGCATGGCCGTCTTCTTTCCCAAACCTCACTGTGGGGAGATGGGTGAAGACAAGT
CAGGCCTTGTTAAAGTTAGTTTCAGAACAATTACTCATGCCTTCCTTTCTCATCCCTAAACATTGGTG
GGGGAGCTACACAATGTACTTTTTCTTTCTAGAGGAAGTATCTATTCATGTGAAAATCTGAAAAATA
TAACAAAGTATGTGAAGATAAAAACCCCTTGCTATTTCA
ACGCGTAAGCGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_017700.2
Summary:	May act as a guanine-nucleotide releasing factor.[UniProtKB/Swiss-Prot Function]
Locus ID:	54848
MW:	25.1