

Product datasheet for **SC208466**

SEC22C (NM_004206) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SEC22C (NM_004206) Human 3' UTR Clone
Symbol:	SEC22C
Synonyms:	SEC22L3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_004206
Insert Size:	677 bp
Insert Sequence:	>SC208466 3'UTR clone of NM_004206 The sequence shown below is from the reference sequence of NM_004206. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
TGGTTCTGCTGGTTGGACAAACCTCGTAGCCAGCCACCCCTGACCCAAATGAGGAGAGCTCTGATTC
TCCCATCCGGGAGCAGTGATGTCAAACCTCTGCTGCTGGGAAATCTCATCAGCAGGGAGCCTGTGGAA
AAGGGCATGTCACTGAAATCTGGGAATGGCTGGATTGCGAAACATCTGCCCATGTGTATTGATGGCAGA
GCTGTTGCCACAAGCGCCTTTTATTTAGGGTAAATTAACAAATCCATTCTATTCCTCTGACCCATGC
TTAGTACATATGACCTTTAACCCCTACATTTATATGATTCTGGGGTTGCTTCAGAAAGTGTATTTTCATG
AATCATTCATATGATTTGATCCCCCAGGATTCTATTTTGTTAATGGGCTTTTCTACTAAAGCATAAAA
ATACTGAGGCTGATTTAGTCAGGGCAAAACCATTTACTTTACATATTCGTTTTCAACTTGTCTGTTCA
TGTTACACAAGCTTCTTACGGTTTTCTTGTAACAATAAATATTTTGAGTAAATAATGGGTACATTTTAA
CAAACCTCAGTAGTACAACCTAACTTGATAAAAGTGTGTAATAATGTATAGCCATTTATATCCTATGT
ATAAATTAATGAGGTGGCTTCAGAAATGGCAGAATAAATCTAAAGTGTATTATTA
ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_004206.4
Summary:	This gene encodes a member of the SEC22 family of vesicle trafficking proteins. The encoded protein is localized to the endoplasmic reticulum and may play a role in the early stages of ER-Golgi protein trafficking. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jan 2011]
Locus ID:	9117
MW:	26.5