

Product datasheet for **SC208467**

LENG4 (MBOAT7) (NM_001146056) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: LENG4 (MBOAT7) (NM_001146056) Human 3' UTR Clone
Symbol: LENG4
Synonyms: BB1; hMBOA-7; LENG4; LPIAT; LPLAT; LRC4; MBOA7; MRT57; OACT7
Mammalian Cell Selection: Neomycin
Vector: pMirTarget (PS100062)
ACCN: NM_001146056
Insert Size: 659 bp
Insert Sequence: >SC208467 3'UTR clone of NM_001146056
 The sequence shown below is from the reference sequence of NM_001146056. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CTTGCCCCGGAGAAGCTCCGGGAGGAGTAAGCTGTACGACGCTCCCTCTGCCAGCTGGTCCCGGGAAT
TCTGTGAACCAGGCTGCTGTCTCCTCCCCAGAAAGAGTCTTACCTTGGAGAGGGTCTGGAGAGAATT
TCCTCTTCCCCAGCTAAATACCCTGCCTGCAACTGAAGCAGACCCGGGGTGTCTCCCTGCCCTCTGC
CCAGAGGCCACCTCCACTCCTACAAAATCAAAGTATTGTCCAGACAAGAGTCACTGGCCCTGCTCCAG
TTCTGGGTATCCAGAGAGCACTGCACTTCCCCAAAACGGAAGGGGCCCTGGGCAGTGGGTTTTGGGC
AAATCCCTTTTCTTGCATCCACAATGTGGGGTCGGAGCTTGGGGCAGGTCTGGGAGTGGGAAGCCT
TTCTTGTGTCTTTGCTCCACTTTTAGCTCATCGCACCAATATTGCAGACTTGGAAAGGAAGCATAAG
TTCCCATTTTCAAAAGGGGAACTGAGGTGCGGGTGCAGGGCTGGGGACGGCCGTCCCATGGCTTC
CATCTGAGCCACCTCGGGACCCAGCACTCCTGGCGCCCTCTTCTCATCGCTTGGCCTATGACAGGTCA
CCGTGTGTAATCTTTCCAATAAAGTGTGCACAAAG
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI
OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001146056.3
Summary:	This gene encodes a member of the membrane-bound O-acyltransferases family of integral membrane proteins that have acyltransferase activity. The encoded protein is a lysophosphatidylinositol acyltransferase that has specificity for arachidonoyl-CoA as an acyl donor. This protein is involved in the reacylation of phospholipids as part of the phospholipid remodeling pathway known as the Land cycle. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2009]
Locus ID:	79143
MW:	23.5