

## Product datasheet for SC208498

### PRR13 (NM\_018457) Human 3' UTR Clone

#### Product data:

|                           |                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | 3' UTR Clones                                                                                                                                                                                                                                      |
| Product Name:             | PRR13 (NM_018457) Human 3' UTR Clone                                                                                                                                                                                                               |
| Symbol:                   | PRR13                                                                                                                                                                                                                                              |
| Synonyms:                 | TXR1                                                                                                                                                                                                                                               |
| Mammalian Cell Selection: | Neomycin                                                                                                                                                                                                                                           |
| Vector:                   | pMirTarget (PS100062)                                                                                                                                                                                                                              |
| ACCN:                     | NM_018457                                                                                                                                                                                                                                          |
| Insert Size:              | 609 bp                                                                                                                                                                                                                                             |
| Insert Sequence:          | <p>&gt;SC208498 3'UTR clone of NM_018457</p> <p>The sequence shown below is from the reference sequence of NM_018457. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> |

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
TCCTCCTCTTCCAGCAGTGATTCTGACTGAATACAGGCCCTGGACCCTTCCCTCAAGTCTCACCAGTTT
TGCTCTCCCATCAAGCTTCAGATGCCATGTTGTACTGGGGGAATGTAGCCCTTGTGCTCCCCACCCCT
ACCTCCACCTGAGCCTCACCTGCTGTTGAGCCCTGAGTGGCTAGGGGAAATGGGAAGAGGATTGCCAT
GGCCTGGCCATCTTGTTGCTGCTTGTTAGATCATATAGCTAATGAATTAGGCAGGGGAGCTATTTTTT
GAAGATGATGAATAAATGTTGAAGACAAGTTTGAGATCTGTAAATGTGATTTTTTACTTCCACTTAT
AATACTTGTGATTGGGGAGGTTTGTGGAATTCATTATGATGAAAAACCTATCTTTTTTGAATGTTG
GCATACTTGGGAATTTAGTGGCAAATACATTCACCCAGCAGGCCTTTTGTGTTGCACTAACTGCAAG
GTTGCTGGGAAGTAGAGTCCATTTGTTGATGAGCTTTGACTGCGGTTTTGGAACCTTACCTCTCTCTC
TTAGCCCAATATGCTGTCTTGGGTCCTATTCAAATAAAGTTATTTCTCCTGGTCTCA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

|             |                                                                                                                                                                                                |
|-------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Components: | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials. |
| RefSeq:     | <a href="#">NM_018457.4</a>                                                                                                                                                                    |
| Summary:    | Negatively regulates TSP1 expression at the level of transcription. This down-regulation was shown to reduce taxane-induced apoptosis.[UniProtKB/Swiss-Prot Function]                          |
| Locus ID:   | 54458                                                                                                                                                                                          |
| MW:         | 22.2                                                                                                                                                                                           |