

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_152994.5
Summary:	Putative regulator of protein phosphatase 1 (PP1) activity (PubMed:19389623). May play a role in the endosomal sorting process or in endosome maturation pathway (PubMed:30520571).[UniProtKB/Swiss-Prot Function]
Locus ID:	129285
MW:	25.6