

## Product datasheet for **SC208637**

### WDR20 (NM\_144574) Human 3' UTR Clone

#### Product data:

|                           |                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | 3' UTR Clones                                                                                                                                                                                                                                      |
| Product Name:             | WDR20 (NM_144574) Human 3' UTR Clone                                                                                                                                                                                                               |
| Symbol:                   | WDR20                                                                                                                                                                                                                                              |
| Synonyms:                 | Bun107; DMR                                                                                                                                                                                                                                        |
| Mammalian Cell Selection: | Neomycin                                                                                                                                                                                                                                           |
| Vector:                   | pMirTarget (PS100062)                                                                                                                                                                                                                              |
| ACCN:                     | NM_144574                                                                                                                                                                                                                                          |
| Insert Size:              | 677 bp                                                                                                                                                                                                                                             |
| Insert Sequence:          | <p>&gt;SC208637 3'UTR clone of NM_144574</p> <p>The sequence shown below is from the reference sequence of NM_144574. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> |

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGCATCGCC
CCTGGTAAAGTGGTAAGTTTAAATCCTTATGCTGCACCAGATCTAGAACTTGAATAGGTAGTGACTTT
TTTCTTTTTCGTGGGAGGGGTGGGGTGTACAATGAATGTGAATGACACTTCTTATTCTTAATGTAAATC
TCAATGCATCAGAGCCATAATTTTGGATACTGCATGCCATGTAATTCTGAATCATTGATAATTTACCT
TAGAGCATTTAAAAAATAATCAAATAATTGCCAGCCAAGTCAGTCATCCTCCTGGGAGTATATAG
AGTCCCAAGGTTAGCGCTCCTGTATTAGACTATTTCAATTTTAGGAAAATCATGACCATGTGGGGAAC
AATGACTTTAAAAATGCTGAAATTTAAATTTATGCTTTAACTGGAATATTTTTTGCTTAACTACTCAATT
AGAATATTGTACACCTGATCAATGTGTGTTTCAGCACAGATGGCCATGAATTGTCATTTATAGTCCAATT
TTTTATCTTAATCATAAAATGTTTAGGAATCTATGAAATTTAACTTTAGGAACAAACGTTTAGCAGGG
TTGATTGATATTATTTTACATTGTTCTGGCAATCCACAGAAAGAGAAGAGCCTTAATTTTTAAACCC
ATTTTAGTCATTTTATGACAATTAAGTTGTTTAATAAACATCTTTTTTCAAAGAA
ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

|                    |                                                                                                                                                                                                                                                                                            |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Restriction Sites: | Sgfl-MluI                                                                                                                                                                                                                                                                                  |
| OTI Disclaimer:    | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs). |


[View online »](#)

|             |                                                                                                                                                                                                                                                                           |
|-------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Components: | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.                                                                            |
| RefSeq:     | <u><a href="#">NM_144574.4</a></u>                                                                                                                                                                                                                                        |
| Summary:    | This gene encodes a WD repeat-containing protein that functions to preserve and regulate the activity of the USP12-UAF1 deubiquitinating enzyme complex. Multiple alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jun 2011] |
| Locus ID:   | 91833                                                                                                                                                                                                                                                                     |
| MW:         | 26.4                                                                                                                                                                                                                                                                      |