

Product datasheet for **SC208690**

zinc finger protein 655 (ZNF655) (NM_001085367) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	zinc finger protein 655
Synonyms:	VIK; VIK-1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_001085367
Insert Size:	699 bp
Insert Sequence:	>SC208690 3'UTR clone of NM_001085367 The sequence shown below is from the reference sequence of NM_001085367. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC TGCTCGGAAACTCCTTTCCACGATTGAGATATTAATAATTGACTGATGGAATAGAAGCTCCCCAGGA TGCCACCACTGTGTAATCGCAGCTCCTCAAATTACCTCTGTTTAATTTCAAATGTTAGGGTCCAAGG AAGCCCTCTGTTGCAACCAGATATGTTTTGAACCCAGTTCATTAGAAACCATGGTTGGTGGTCATCAT CTACTTGTATTGTGAAAAACCAGAAATTCAAATTCAGCTCTTCAGTTATGGTCTTTTCATATTATTTT AGAATGTGAAACTTATGCTGTGTGCCCTCCTTTATGTGGCGAACACAACCTGCCCCGTATGGTGCCT AACTCTCCTGTGAACCCTGGGCTTCCTCATGTTAACATTCTGGTTTCCTCCTTGAAGTCCCTG TGGCCCTCTCAAGTCTACCTTCTACCCTTTATTTCTTGTGGTCTGTACAGATTGATGCAATTAT ATTGCCTTTCTGATTAAGTGTTCCTTACTCTGTCTTGTATGGAGGGGAGAAGGGAGGCAAGAAGTT AAATTTATTTCAAGGCATCAAATTGAATATCTACTTAATAAATTTGTAAGTGTAAAGTGACATACATA TAATTGTAAGTACTGTTAATTTTATCTAATAATTACAGAGTTTCAACAGTAAAAAAAAAAAAAAAA AAAAAAAA ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001085367.1</u>
Summary:	This gene encodes a zinc finger protein. The zinc finger proteins are involved in DNA binding and protein-protein interactions. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
Locus ID:	79027
MW:	26.9