

Product datasheet for **SC208701**

zinc finger protein 655 (ZNF655) (NM_024061) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	zinc finger protein 655
Synonyms:	VIK; VIK-1
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_024061
Insert Size:	676 bp

Insert Sequence: >SC208701 3'UTR clone of NM_024061
The sequence shown below is from the reference sequence of NM_024061. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
TGCTCGGAAACTCCTTTCCACGATTGAGATATTAATTTGACTGATGGAATAGAAGCTCCCCAGGA
TGCCACCACTGTGTAATCGCAGCTCCTCAAATTACCTCTGTTTAATTTCAAATGTTAGGGTCCAAGG
AAGCCCTCTGTTGCAACCAGATATGTTTGAACCCAGTTCATTAGAAACCATGGTTGGTGGTCATCAT
CTACTTGATTGTGAAAAACAGAAATTCAAATTCAGCTCTTCAGTTATGGTCTTTTCATATTATTTT
AGAATGTGAAACTTATGCTGTGTGCCCTCCTTTTATGTGGCGAACACAACCTGCCCCGTATGGTGCCT
AACTCTCCTGTGAACCCTGGGCTTCCTCATGTTAACATTCTGGTTTCCTCCTTGAAGTCCCTG
TGGCCCTCTTCAAGTCTACCTTCTACCCTTTATTCTTTGTGGTCTGTACAGATTGATGCAATTAT
ATTGCCTTTCCTGATTAAGTGTTCCTTACTCTGTCTTGATGGAGGGGAGAAGGGAGGCAAGAAGTT
AAATTTATTTCAAGGCATCAAATTGAATATCTACTTAATAAATTTGTAAGTGTAAAGTGACATACATA
TAATTGTAAGTACTGTTAATTTTATCTAATAATTACAGAGTTTCAACAGTTAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_024061.4</u>
Summary:	This gene encodes a zinc finger protein. The zinc finger proteins are involved in DNA binding and protein-protein interactions. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
Locus ID:	79027
MW:	25.9