

Product datasheet for **SC208717**

ZNF83 (NM_001105549) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: ZNF83 (NM_001105549) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: ZNF83
Synonyms: HPF1; ZNF816B
ACCN: NM_001105549
Insert Size: 667 bp
Insert Sequence: >SC208717 3'UTR clone of NM_001105549
The sequence shown below is from the reference sequence of NM_001105549. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GCCGAAAAGAAATCTAACACATGTAATTAATGTGGCAGAGTGTTCAAGTTAGCATTAAAGCCTTGTAAGA
CATACAATAATTTATACTGGAGAAAACTTTGCAAGTATAATGAATGTAGCAGAGCCTTTAGTTTTGT
TCAAGGCTTAATAACCGTTAGCTAGACCATAGAGGACAGAACTTTACTAATGTAAGTGAATGTGGCAAG
GTCTTAAGGTAAAATCTGAGACCAGGATTCTTCAAAGAATTCTTGCTGGTGAGAAACCTAACAAATGTA
ATGAATGTGGCAAGTCTTCTGGCACAATTCTCACATTGTACAATATTGCAAAAATTCATGCTTGAGAG
AAACAAAAACACTGAGAGTGGGAAACCATTATGACTTCAAACATTCATCAACATCAGAGAATCCATACT
AAAGAGCATTTATAATAATTATATGTGATAGAGATTTCCGCGAGGCCAAAGTCTCACTAGGCATCAAAA
ACTTTTTTGATGAAACCATACAAATGTAACTGATGCTTAAGCTTTTACCCAGGCATCAAAACCGGAA
CATCACAGGGTTTATACTGGAGAGTAACACAAAGATAATGTAATAAGCCTTTTCAGTGAATATTCA
TGATTTTGTGCTGAGAGATCCACTCAATAAAAACCGGCAATGTA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



[View online >](#)

RefSeq: [NM_001105549.2](#)
Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID: 55769
MW: 25.9