

Product datasheet for **SC208729**

CCDC92 (NM_025140) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	CCDC92
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_025140
Insert Size:	1462 bp



Insert Sequence:

>SC208729 3'UTR clone of NM_025140

The sequence shown below is from the reference sequence of NM_025140. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAGCACTCAGGGACGGACAGAACTGTGTGAAGCCCGCCGTGCCCCACCCCGCGCTGTCCATGCACTGTG
AGCACCACCTGGGAAATCTCAGCCACACCTTTTCTGTTAATCCCATGCATGCCAAACACTTTTCACACC
TACCGACCCATTCTCCTTCTGCTTCTCTTGCCCTCTTCTTACACCAAAATATGATCGTGTCCCTGCCG
CAGAATATGATTTTCTAATTGCTGTGGCCAAACGCCTGTGTGCCGAATCGTTGCTTCTGATCCCGCT
CCGTGTAACCTAAGTGCCTGCAGGCAAGCCCAAGGCCACGGCTGCGTCACTACTGATGTTTACGATGC
CACACAGTCACACACCTAATTCATTCTCAAGTCGGAGCAACACATACCAACCTTGACCTTATCCTCAAG
CTCCAGGGCAGCCTGGCCGAGCAGCCCTGCTCCCTCTGGAGACCCTGTGACCTCCCGAGCTCCTCC
TGGAGACCCCTGTACCTCCTGACCAACCTTTCCAGGGCGGCACCGATCACCGAGCAGCCGTGCGTGT
ATCTCAAGGAATAAATAAGATGACGCTACTCCTCATAGCACCAACCTGAATGTGTGTTTATATTTT
TTTGTTAGTTTTATCCAAATGTTTAAAGATCCCAACAACTTTATTTTCTAAACCTGCCTTCTCTGT
TTTGTGTCTTTTATTAACAAACACCTGAAGGAGGGCTAGGAGACTGCAGGAGGGAGAGCAATCAGAG
TCCAGAAGCGACTCCGGCCTCCTCCCGCCAGAGCTAGTCCTTAGAATGTGCTGCGCTCCTGAGCTT
CTGTGGGAAGCTGCTTTTGCCTCTGATGGCCCTGCCCACTGTGCTGCCCGTGTAGCCCTCTCCTCC
CTCATGGAGTGAAGCCACCCCTCACCTCCAGGTGGCCAGCCAGGGCTCCAGGCCAGCCTCATGGTGCAC
AGAGGGAAGGGGTGCCAGTGACCACTGGTCTGGAGTTCTGTGCTTACCTGACCCACACAGACC
CAGGCAGGCTGGGGGTAGAAGCTCATGCAGGAATAACACTAAGCTCAACAACTTCACAGTCGCCGTA
GCAGGTGAGAGCCAAAGAAGCTAGGGCCAGAGAGTTCTGGAAGGGGAGAGTGGCTCCAGGTGACTGG
TCATTCTGCTCTGTGCCGACATCATAGCCCTTCTTGCCAAAAGGAGCATCTGCTGAAACAGACCT
GGATCATCGAAAACACACCTTTGGAATTACATTAGGATTCTCACTTGTGCCAATGCTATGTATGCCAG
GCATTTTGCCTAATCAAGTGTCTCTGAGAAGGAAAAGAACTCCACACATCTTAAATACTCTC
CAATATTACTTTTTGATTATATAAATTCCTTTTAAAAAGTAGTATTAACAAACAAACCATTCAGGTTA
ATTTGTATGGAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites:

SgfI-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

[NM_025140.3](#)

Locus ID:

80212

MW:

54.2