

Product datasheet for **SC208732**

C10orf30 (BEND7) (NM_001100912) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	C10orf30
Synonyms:	C10orf30
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_001100912
Insert Size:	685 bp

Insert Sequence: >SC208732 3'UTR clone of NM_001100912
The sequence shown below is from the reference sequence of NM_001100912. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GGCATTGCTCTACCAACAGTGGTCTGATGAAGGACTGCAAAACACAGCTCTTGTGTTTTTGCACC
ATGGGTTTTAAATCCGTCCTCTGCTCCCAAGGAGCTACTACCAAGGGAAGGAGAAGATTAATT
GGTGACCTTTCGCGGTGTTAAATAGGATGAGTATGTTGCTGGTCAAGGCAGAAACGCTAACTCCAAC
TGTGGTCTGTAAAGAGATCTGCTGGGAACCTTCAGAAGCACTGCCTTCCAAGGCGGGGCCACAGCTTGG
GGTGTTCAAGCACTTGAGGACGGGTGGAGCTTGTTCACCCAAAGATGTGTGCGATGAACGGATCTACAC
ACTGGAATACGTATCAAACTGCAACAGCAAAATCATATTTGTAAATAAATTTGTGTTACTTACT
GTGACTACTGTTCACTCTTTATTCAGAAACCTTTTATAGGCTTTCTAGCAAAAAAAAAATCCATATC
TATGGATAATGTCTCAATCTTCTGTATATATGTTTTATTAATATGTAATATTTAGATTTTTTAAAT
ACTATGTTCTTTAAAAAAGAATTCAATGCAAGCCTGGTTGTGAAATGGTGTATAACATTGTGTTGTGG
TATCAATCCCAAGATGCTCTTTATGTCTGATCTGGAAGAATAACAATAATTACTTTAATTGAA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001100912.2</u>
Locus ID:	222389
MW:	26