

Product datasheet for **SC208747**

TPD52L1 (NM_001003397) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	TPD52L1
Synonyms:	D53; TPD53
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_001003397
Insert Size:	1594 bp



Insert Sequence:

>SC208747 3'UTR clone of NM_001003397

The sequence shown below is from the reference sequence of NM_001003397. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AAGAAGTTCGAGACATGAGACGAAAGTAGGCGGTACGAACCTAATGGAGGCAGTTTTGAGGAGGTCC
TCAGCTCCACGGCCCATGCCAGTGGCCAGAGCTTGGCAGGAGGCTCCCGCGGACCAAGGAGGAGGAGC
TGCAGTGCTAAGTCCAGCCAGCGTGACAGTGCATCCAGAAACCGGCCACTACCCAGCCCATCTCTGCCT
GTGCTTATCCAGATAAGAAGACCAAAATCCCGCTGGGAAAAACCCAGGCCTTGACATTGTTATTCAAAT
GGCCCTCCAGAAAGTTAATGATTTCCATTTGATTTGTGTTGATGATGGACCACTTGACCATCACAT
TTCAGTATTCATAGATGACTGTCACATTTTAAATGTTCCCACTTGAGCAGGTACACAACCTGGTCATAA
TTCCTGTCTGTGAATTCGATGTATATTTTCCAAACATGTAGCTATTGTTTGTCTTTGATTTTGTCTTG
GCCTCCTTTATGATGTGCATGTCCTTGAAGGCTGAATGAACAGTCCCTTTCAGTTCAGCAGATCAACAG
GATGGAGCTCTTCATGACTGTCTCCAGCAATAGGATGATTTACTATAAAATTCATCCAACCTACTTGTGA
TCTCTCTCACCTACATCAATTATGTATGTTAATTTAGCAATTAAGAAGATTGATTTTAATGACTTTGA
ATTCTTAATTTCTTTGTCTTAAAGTTGCTAGTTATGATTTTACAGATGCAATTTTAAATCAACTTTTA
GCCAGGTGCGGCGGCTCACACCTGTAATCCCACTATTTTGGATGCCAAGGTGAGAGGATTGCTTGAGG
CCAGGAGTTAAAGATCAGCCTGGGCAACACAGACCCTGTCTCTACAAAAAAGAAAAAATTAGCCAGA
CATAGTGTTGCTTGCCTGTAGTCCAGCTACTCTAGAGGCTGAGGCAGGAGGATTGCATGAGCCTAGGA
GTTTCAAACTGCAGTGAGCTATGATTGCACCACTGCACTACTCCAGCCTGGGTGACAGAGTGGGACACT
GTCTCCAAAAATAGTAATAAAGTAGTCAACTTTTACTGCTAATTTGGTGAACATGAGAGAGGATATG
AAAATAAATATTACCTCAGCTATCCTAGGATGTTAAAAATACTCCAATTTTAAATTTCTCTCCAATCT
ACATACAGTAGTAGTTAGTCAGATAAAGGATATCCAAAAAGAGATAGCTAGAAAAATGGGAGAAGCAGA
GTTCTGCAACCCCTTTCAGTTTGTAATTTGCATGTATGAAAAATACTGGTATTATCAATCCACT
CAGATTTCTGCACTAACTTTTATCTTATATCATATGTATCTCTTTCTTTTCTAAATGGGAACATA
TATTTGTTATTAGGTGGCAGAGATATAGCCTTAAGATATATTTGTAATAATGCACACTGAATAGACATCC
AACCTAAAAAATCACTATTTAAAAAGCCCATATAATATATACATATTTGTTAGCATGCTAATTGTTT
ATGTTTTGTGTTTATTAATAGAAGTGATATATGACATTTTGAAGTAAAGCACATCTGAAAAATTCT
ACTCAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCACCGCCGCTCTATGAAAGG
```

Restriction Sites:

Sgfl-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

NM_001003397.3

Summary:	This gene encodes a member of a family of proteins that contain coiled-coil domains and may form hetero- or homomers. The encoded protein is involved in cell proliferation and calcium signaling. It also interacts with the mitogen-activated protein kinase kinase kinase 5 (MAP3K5/ASK1) and positively regulates MAP3K5-induced apoptosis. Multiple alternatively spliced transcript variants have been observed. [provided by RefSeq, Jan 2016]
Locus ID:	7164
MW:	61.1