

Product datasheet for **SC208765**

THUMPD2 (NM_025264) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: THUMPD2

Synonyms: C2orf8

Mammalian Cell Selection: Neomycin

Vector: pMirTarget (PSI00062)

ACCN: NM_025264

Insert Size: 704 bp

Insert Sequence: >SC208765 3'UTR clone of NM_025264
The sequence shown below is from the reference sequence of NM_025264. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TATAAGAAGTCGCACTCTTCTGGACTGTAGCAGGCTTGCTGCCATCAGCCGGGTTTCAGGCCCTTGTAAAG
TTAGCTGTACAGCAGAAGTTTGTCTCTAGGATTCATATAGCACCCACAGAGCTCCAAGTAACCAAA
TTCCCCAAAGACAGGAGGTGTGGCTGATGGGTACACATGCTTGGTAAGAGTTTCTGCTTTATATTCTA
AAGTACTTAATGTTTTTGAGAAAAAAGTTATCGAATATTTTTTAGACTTAAGTTTTATCTTCTGTCAGA
ATGTGAAAAGTTCTGCCAGTTTCTGAGCAAAAACCTCAATTTCAAATAAAGGTTTTTCAGGGGTTTAC
TTAAGTGTATCACAATCAACACCAGACGGGATATTTAACTATCAGATTTTCGAAATGTAAAAACCTT
TATGTGTAACCTCTTTGATATTAACCTTGATTTTGGCCTGTAAATTTCTTCTCCATTGGAATGAAA
ATAATTTGTATTTGCAATACATTTTATATGTTTATTATGTAACCAACAGTGAGAATAATACTGAATT
AAGCCAGTATCCTTCAAACATGATGGAATAAATTTGCTTTTTTCTAAAAGCATTATACACTGT
CACCTCGCCCTTATTTATAGAAAGTCAGTTTGTGCTTATGTTCACTACCTGAGAAAAAGAAACAGT
AAGTTTGAGGAGAA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_025264.5</u>
Locus ID:	80745
MW:	27.3