

Product datasheet for **SC208772**

GIPC (GIPC1) (NM_202470) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	GIPC (GIPC1) (NM_202470) Human 3' UTR Clone
Symbol:	GIPC
Synonyms:	C19orf3; GIPC; GLUT1CBP; Hs.6454; IIP-1; NIP; OPDM2; RGS19IP1; SEMCAP; SYNECTIIN; SYNECTIN; TIP-2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_202470
Insert Size:	684 bp
Insert Sequence:	>SC208772 3'UTR clone of NM_202470 The sequence shown below is from the reference sequence of NM_202470. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
ATTGGGGACGCCAAGGTCGGCCGCTACTAGGACTGCCCCGGACCCTGCGATGATGACCCGGGCGCAAC
CTGGTGGGGGCCCCAGCAGGGACACTGACGTCAGGACCCGAGCCTCCAGCCTGAGCCTAGCTCAGCAG
CCCAAGGACGATGGTGAGGGGAGGTGGGGCCAGGCCCCCTGCCCGCTCCAATCGGTACCATCCCCTCC
CTGGTTCCAGTCTGGCCGGGTCCCCGGCCCCCTGTGCCCTGTCCCCACCTACCTCAGCTGGGTCA
GGCACAGGGAGGGGAGGGATCAGCCAAATTGGGCGGCCACCCCGCTCCACCACCTTCCACCATCAGC
TGCCAAACTGGTCCCTCTGTCTCCCTGGGGCCTTGGGTTCTGTTGGGGGTCATGACCTTCTAGTTTC
CTGACGCAGGGAATACAGGGGAGAGGGTTGTCCTTCCCCCAGCAATGCAATAATGCCCTCACCCCTC
CTGAGAGGAGCCCCCTCCCTGTGGAGCCTGTTACCTCCGATTTGACACGAGTCTGCTGTGAACCCCGC
AACCTCCTCCCCACCTCCCATCTCTCCTCCAGGCCATCCCTGGCCAGAGCAGGAGGGAGGGAGGGA
CGATGGCGGTGGGTTTTGTATCTGAATTTGCTGTCTTGAACATAAAGAATCTATCTGCTGTT
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_202470.3
Summary:	GIPC1 is a scaffolding protein that regulates cell surface receptor expression and trafficking (Lee et al., 2008 [PubMed 18775991]).[supplied by OMIM, Apr 2009]
Locus ID:	10755
MW:	24