

Product datasheet for **SC208810**

SFRS4 (SRSF4) (NM_005626) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: SFRS4

Synonyms: SFRS4; SRP75

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_005626

Insert Size: 697 bp

Insert Sequence: >SC208810 3'UTR clone of NM_005626
The sequence shown below is from the reference sequence of NM_005626. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TCTAGATCTAGGTCCCACTCAAGGTCCAACTGGCTATGGCCACAGCTGGAACCTACCCGAGAAGTCTTT
TGTACATGTTTGGTAGCCGTAGCACAAAGTGATTGGAGTAGAACATGTCAGTGTACATTTTAACTC
CCCTAATGGTGTGTCTATAATTGTTAAATCTAAGTGCTTCTCTCAGTAAAGCCTCCTGGCACCAGGCC
TTCCTGCTCGACTGAAAAAATTTTCTTTGAAAAATCCCCTTTTACTCATGGCCACAGTAGAATATC
CAAAACGCCTTGGCTTTTCAAGCCTGGCCTTTCTACAGGGAGCTCAGTAACCTGGACGGCTCTAAGGCT
GGAATGACCACATAGGTAGGTATGGTGAGTTCAACCATTTTGGCTCTTGAATTGATGCCCTTCGATGTA
TGCCATTTAGTGAAAGTGCTAAGTCTTAAGTTTCTACCACTTTGGTTTCATATTTTGGACTTAACAA
AGTTGTGAATAGCACAGTCGAGGAAAATTGATACCTGCAGTAACCCATAGGAAATAAACTGTAGAGTTC
CATATTCTGGTATTGTGATTATATTGTTTTATATTAAGAAAGAAAGAAATTTTTTTAATTTTA
TTTTTCCCGTCTTGCAAAGTATAGTGACCCCTGTTTCCATTAATTTGAATAAAGACTATTTTTGCTT
GACAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAAGAAATGACCGACCAAGCGACCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_005626.5</u>
Summary:	This gene encodes a member of the arginine/serine-rich splicing factor family. The encoded protein likely functions in mRNA processing. [provided by RefSeq, Feb 2009]
Locus ID:	6429
MW:	26.3