

## Product datasheet for **SC208832**

### **TTC13 (NM\_001122835) Human 3' UTR Clone**

#### Product data:

|                  |                                                                                                                                                                                                                                                          |
|------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:    | 3' UTR Clones                                                                                                                                                                                                                                            |
| Symbol:          | TTC13                                                                                                                                                                                                                                                    |
| Mammalian Cell   | Neomycin                                                                                                                                                                                                                                                 |
| Selection:       |                                                                                                                                                                                                                                                          |
| Vector:          | pMirTarget (PS100062)                                                                                                                                                                                                                                    |
| ACCN:            | NM_001122835                                                                                                                                                                                                                                             |
| Insert Size:     | 692 bp                                                                                                                                                                                                                                                   |
| Insert Sequence: | <p>&gt;SC208832 3'UTR clone of NM_001122835</p> <p>The sequence shown below is from the reference sequence of NM_001122835. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> |

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
TCTTCTCCACGTTGTCTTAAGAACTCAGTTCTGCTGCTGTATTTATACAAGTAAAGGGCCGGACCTC
TTGCTTCTTAAGTTATTTTTTAAACATGGAATTATAAGAAATTAGGTCCTCTATTACTTTTGATACC
AATTTTATAGGAATTGAATACTTGAATCATTCTTTACTTTTCTGAACCATAATGAAAATCATGAA
AGAGGCTTTTCGGGGGAAAAAGCTACTTGACTTGAGGGGAAAGCGTATGTATTAGCCGGCAGCCTGTTG
GTTCCCACTTATAAAGTCATATAACTCCTTCCTTACAGTAAAAAAATTTTCTCAACAATAAAATAT
AATACTGAAGATGAACAAAAACAAACCAATAGGAGTGTTTGTATTCTTTAGAGATCTGGATATTTT
TCAAAGACTTGCTTTTACGCACATTGTGATAGAGCATCTTTTGTGACTCACGTAACCTCCATGTGTG
GGTCTTTGTTTATTAGAACCATGGAAGAAGTTCTGGGCCATCAGTATGTGTGCAGAGCTCACAG
CCAGGAACACTTATTTGAACAATGCAGACGACTTGCCAGAGTCCCAATGCTGAATTATTGAACCA
TGTTTTTTTTCTTTGTTGAAAAAATAAACAGGGTGATTCTATAAATCACCAGTTTTTGGCATTG
GA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCCTTCTATGAAAGG

```

**Restriction Sites:** SgfI-MluI

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



|                    |                                                                                                                                                                                                |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Components:</b> | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials. |
| <b>Note:</b>       | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                       |
| <b>RefSeq:</b>     | <u><a href="#">NM_001122835.3</a></u>                                                                                                                                                          |
| <b>Locus ID:</b>   | 79573                                                                                                                                                                                          |
| <b>MW:</b>         | 26.7                                                                                                                                                                                           |