

Product datasheet for **SC208838**

ZFP64 (NM_199427) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: ZFP64 (NM_199427) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: ZFP64
Synonyms: ZNF338
ACCN: NM_199427
Insert Size: 576 bp
Insert Sequence: >SC208838 3'UTR clone of NM_199427
The sequence shown below is from the reference sequence of NM_199427. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GCTCCCCTAGAGCCCAGCCAAGACCTCTAGCTCAGCCGAAGAGGGTGGTCAGCTGTCCGAACTTCTGAC
CAGTGCTGTTGCCAGCCCTCCAGAGGTGAGGCCACTGAGAGACAGGCTTTGGCACCACCTCTGGGGC
TGGCCGCTCCAACCTTACAGGCCTTCAGAGGCGATGGCTGGGCGACAGTGACAAAAGCAAAGCAAAGCA
GGGCTGTGGAGACTCCTCGCATTGTCTCTTCCCTCCAAGGATTATCTGAGCAAGTCGACTTGTTCAT
TTCAAAGCGGGGTCTGCCAAGCCCTGCTCTATCCAATGGGGATAGCTTCTACGTAACGGATTCCAATT
TAAAACTGATCCATCTCTGATGGAATTATAACTTAGACGATATAGAAAATAAACTAATCTCTAG
AGTGAAAGGGTTAAAGTGAGGAAGGTGTACGTGAAGGGTCTGGGGTGTCTGGAAAATTCTATTTTCAT
AACCTGGGTGATCGTTCTAAGAGTGTTCCTTATATGTTGTGTGGTTCTTTTACCTGTTGCTTAC
AATAAAAAAGGCAGACTCATGGTA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 µg dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



[View online »](#)

RefSeq: [NM_199427.3](#)

Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Locus ID: 55734

MW: 21.3