

Product datasheet for **SC208849**

NEK5 (NM_199289) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: NEK5 (NM_199289) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: NEK5
ACCN: NM_199289
Insert Size: 681 bp
Insert Sequence: >SC208849 3'UTR clone of NM_199289

The sequence shown below is from the reference sequence of NM_199289. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CCGTACAGCCCAGTGCTCATCTGATGTAGTAGTCTACTTCTCACTATACACCCTATAGATCTTGTATC
AGACACTTTCAAATATGTTGTTTTGATATCTCCCTATGCCTAAAGAGACTTGCCTGCTAATTATCCTGG
ACAAATGACCACTCTAGACATGCACATATAAAGGAGGAAAAATCCAGGATGCCACCTTCCTCAGAGTGA
GGGAAAATGCGTGAGGAAGGACAGCTGACTTGATTACCTTCCAAGGCCAGCTCAACAACCTTTTTTGC
GATTTTCATAGGGTCATTTTTTTTTTTTTAATTTTTATTTTTATTGATCATTCTTGGGTGTTTCTCGCA
GAGGGGGATTTGGCAGGGTCATAGGACAATAGTGGAGGGAAGGTGAGGAGATAAACAAGTGAACAAAGG
TCTCTGGTTTTCTAGGCAGAGGACCCTGTGGCCTTCCGCGAGTGTGTTGTCCTGGGTAAGTACTTGAGATT
AGGGAGTGGTGATGATTCTTAACGAGCATGCTGCCTTCAAGCATCTGTTTAACAAAGCACATCTTGAC
TGCCCTTAATCCATTCTCTGAGTGGACACAGCACATGTTTCAGAGAGCACAGGGTTGGGGTAAGG
TCACCGATTAGCAGGATCCCAAGGCAGAAGAATTTTTCTTAGTACGGAACAAAATGAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-Mlul

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_199289.3](#)



[View online »](#)

Locus ID: 341676

MW: 26