

Product datasheet for **SC208858**

TMEM234 (NM_019118) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	TMEM234
Synonyms:	AASL548; Clorf9l; dJ622L5.7; PRO1105; RP4-622L5
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_019118
Insert Size:	683 bp

Insert Sequence: >SC208858 3'UTR clone of NM_019118
The sequence shown below is from the reference sequence of NM_019118. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AAGACCCAGGGGCAACAGTCTACCTTTAGTGGGCCGAACCACTTCCAGCTCTGCTGCCTCCAGGAA
GCCCCCTGGGCCATGAAGTGCTGGCAGTGAGCGGATGGACCTAGCACTTCCCCTCTCTGGCCTTAGCTTC
CTCCTCTCTTATGGGGATAACAGCTACCTCATGGATCACAATAAGAGAACAAGAGTGAAAGAGTTTTGT
AACCTTCAAGTGCTGTTTCAGCTGCGGGGATTTAGCACAGGAGACTCTACGCTCACCTCAGCAACCTTT
CTGCCCCAGCAGCTCTCTCTGCTAACATCTCAGGCTCCCAGCCAGCCACCATTACTGTGGCCTGAT
CTGGACTATCATGGTGGCAGGTTCCATGGACTGCAGAACTCCAGCTGCATGGAAAGGGCCAGCTGCAGA
CTTTGAGCCAGAAATGCAAACGGGAGGCCTCTGGGACTCAGTCAGAGCGCTTTGGCTGAATGAGGGGTG
GAACCGAGGGAAGAAGGTGCGTCGGAGTGCGAGATGCAGGAAATGAGCTGTCTATTAGCCTTGCCCTGCC
CCACCCATGAGGTAGGCAGAAATCCTCACTGCCAGCCCTCTTAAACAGGTAGAGAGCTGTGAGCCCCA
GCCCCACCTGACTCCAGCACACCTGGCGAGTAGTAGCTGTCAATAAATCTATGGTAAACAGA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_019118.5</u>
Locus ID:	56063
MW:	24.8