

Product datasheet for **SC208861**

ZC3H18 (NM_144604) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZC3H18 (NM_144604) Human 3' UTR Clone
Symbol:	ZC3H18
Synonyms:	NHN1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_144604
Insert Size:	695 bp
Insert Sequence:	>SC208861 3'UTR clone of NM_144604 The sequence shown below is from the reference sequence of NM_144604. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAGCGGGCCAAGATCCCCGGGAAAGCATAGGCCGTGCCCGACCGGACTGGACGCATTTTTATACATAG
GGTAAGCGCAGCCATTTGGATTTTGCAGTTAATGTCTATTTTGGCTGTGATTCTTTTTAAAAAGTAA
AAAAGAAAAAAGTTTCTCAGCTGAAAAGAAGCCACACAGGAAATGACAACGACGCTGAATCCCAGC
CTCCCTCCCAGAGCAGAAGTCCCGCAGGACAGACAGACAGACAGCGTAGTGACCAGCACGGTTCT
CATGTAAATTACAAGCCCCAGCCGCCAGCCCCGCCTTCTCTTCTCCTCCTCCGTCTTCTCCCTGGCC
CTGGTCAGGCCTGTGGAGCCCCAGCTCTGGGTCCCTAGCCGGGTCCAGGCAGCCAGGCTCCCTCCTGA
GCTGAGAAACGGAACTCGGAACCACTGGTGGCACATCCTTCTCCTCCCCGCCCTGATCACCCGCC
CCCGGATCAGAAATATATCTATATTCTCGACTAAAGTCTCATCAGGAAATATTTCTGTCTTTATTTT
AAGCATCAAATTGTTTAGTTGATTTAAAAAGGAAAAATACAGAAAAGACCAAAAAAGGCCAAGGGT
GTTGTTGGGGCTGTCTAATGTGGTGGGTCTTTTTTGGGGGTCTCCTAAAAATAAATATTTTGTAT
AAGCA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_144604.4](#)

Locus ID: 124245

MW: 25.7