

Product datasheet for **SC208876**

CGREF1 (NM_006569) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: CGREF1

Synonyms: CGRI1

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_006569

Insert Size: 710 bp

Insert Sequence: >SC208876 3'UTR clone of NM_006569
The sequence shown below is from the reference sequence of NM_006569. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ATTGTTCAAGTGGAGAATGATGAGATCAGATCTTGAAGATACAGGTACCCACGAAGTCTCAGTGCCA
GAACATAAGCCCTGAAGTGGGCAGGGGAAATGTACGCTGGGACAAGGACCATCTCTGTGCCCCCTGCCT
GGTCCCAGTAGGTATCAGGTCTTTCTGTGCAGCTCAGGGAGACCCTAAGTTAAGGGGCAGATTACCAAT
AAAGAACTGAATGAATTATCCCCCGGCCACCTCTCTACCCGTCCAGCCTGCCAGACCCCTCTCAGAG
GAACGGGGTTGGGGACCGAAAGGACAGGGATGCCGCCTGCCAGTGTCTTCTGGGCCTCACGGTGCTCCG
GCAGCAGAGCGCATGGTGCTAGCCATGGCCGGCTGCAGAGGACCCAGTGAGGAAAGCTCAGTCTATCCC
TGGGCCCCAAACCCTACCGGTTCCCCCTCACCTGGTGTTCAGACACCCCATGCTCTCCTGCAGCTCAG
GGCAGGTGACCCCATCCCCAGTAATATTAATCATCACTAGAACTTTTGTAGAGCCTTGACACATCAGG
CATCATGCTGGGCATTTTATATATGATTTTATCCTCACAATAATTCTGTAGCCAAGCAGAATTGGTTCC
ATTTGACAGATGAAGAAATTGAGGCAGATTGCGTTAAGTGCTGTACCCTAAGGTGATATGCAGCTAATT
AAATGGCAGATTTTAATCCA
ACGCGTAAGCGGCCCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_006569.6</u>
Summary:	Mediates cell-cell adhesion in a calcium-dependent manner (By similarity). Able to inhibit growth in several cell lines.[UniProtKB/Swiss-Prot Function]
Locus ID:	10669
MW:	25.9