

Product datasheet for **SC208885**

PRRT3 (NM_207351) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: PRRT3 (NM_207351) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: PRRT3
ACCN: NM_207351
Insert Size: 715 bp
Insert Sequence: >SC208885 3'UTR clone of NM_207351

The sequence shown below is from the reference sequence of NM_207351. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGCGCCTCCAGTGATACCATCGAGCTTTGAGAGCGGTCTGACGCAGGGCCAGGACCCTGCCCGATGC
CCACATGGCATGGCTACTGGGACAGTTGAGCCTTTAAAAAATGGGCATCCTGGCCAGGCGCGGTGGCT
CAGGCCAGGAATCCCAGCACTTTGGGAGGCCGAGGTGGGTGGAACACGAGGTCAGGAGTTCGAGATCAG
CCTCCTGGTCTCGAACATAGCGAAACCCCTCGTCTCTACTAAAAATACAAAAAATTAGCCAGGTGGCAC
ACGCCTGTAGTCCAGCTACTCGGGAGGCTGAGGCGGGAGAATTGCTTGAACAGGGAGGCAGAGGTTGC
AGTGAGCCGAGACCACGCCATTGCACTCCAGCCTGGGTGACAAGAGTGAGACTCCGTCTCAAAAAAAAA
AAAAAAAAAAAAAAAAAAGGCATCCCAGGCCTCTACCTGACGGCAACCCCGGGGGACAGCCA
GGTTTGAACACACCCACCCCATGCCGTTTTTTGCTTTAATATTTCTATTTTTTTTTCAGCCGGCAT
TTTGCACAGAGGCCGTGTCTTATGCGGAATACAGGGTGGGTGTGCATGGATTGGGCATGAAGAACAGG
GAGCAGGTGCCCCAGACTCCCCAACGTGGAGTGATGAGCCTCGCTGTCTAGATTGTCCTCTTTGTGCC
AAAGAAAATAAACCTTAGGACTTGG
ACGCGTAAAGCGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.



[View online >](#)

RefSeq: NM_207351.5

Locus ID: 285368

MW: 26.6