

Product datasheet for **SC208901**

TM6SF1 (NM_023003) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: TM6SF1

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PS100062)

ACCN: NM_023003

Insert Size: 766 bp

Insert Sequence: >SC208901 3'UTR clone of NM_023003
The sequence shown below is from the reference sequence of NM_023003. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAAACAAAGGCAGAAGAAAAAGTGAATAAATAATTACTTCATGTTCTCTCTTCTAAATTACTAACT
TTTGTATACTGGTACTGATATTTTGTCCATTTCACCTCTCTCTCATACGTGAGTACTTAAGAATATG
TACATTCCTGCTCTGCACTGTATGTGTGAGCTATATGGTATTGTGTAAATTTTTTTGAAGGAAAATGG
AAATTCCTGAGAAACAGTTTGTGTTAAAGAAATATATCAAAATCATTGTGAATAAACTTGATCATCCA
TCTCAATATTGTTTGACATATAAAATAATTATAAGTGTAATAAATTACAATTTAGTGCCAACAGTAGTG
AGCATGAAATGAACTATTCAAAAGAGAATATGGCCTGTGCATATTAAAAAATCAAAACAGTGAATGC
AGACTGGAGGAGTAACCTTTGCAAATAAGATGAATATGCTTCATTATAACTCAATATAAAAGGCAAA
TCATCAGAATATTTTAAATGTTGTTTGAATAATGTTTCCCAAGGAAAGTTTATTATTTGCTGCTGTT
TCAAGAAAATTACTTTTACTAAATTTTTTGTGTGAATTTAAACAGCTAAATAGGGATCAGTAACCTTA
TCTCTATCCTTAATGAACATTTGTTTTATTGGTGGCTGGAATATTCTATTGTATTCTGTGTATATT
TTTAATAAAATTATTTTGGCCTCTTATAAGACAAATCTTATTAATTTGAAACATTTTCATATATACA
CATAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_023003.5
Summary:	May function as sterol isomerase.[UniProtKB/Swiss-Prot Function]
Locus ID:	53346
MW:	30