

Product datasheet for **SC208902**

TM6SF1 (NM_001144903) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	TM6SF1
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_001144903
Insert Size:	766 bp
Insert Sequence:	>SC208902 3'UTR clone of NM_001144903 The sequence shown below is from the reference sequence of NM_001144903. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAAACAAAGGCAGAAGAAAAAGTGAATAAATAATTACTTCATGTTCTCTCTTCTAAATTACTAACT
TTTGTATACTGGTACTGATATTTGTCCATTTCACCTCTCTCTCATACGTGAGTACTTAAGAATATG
TACATTCCTGCTCTGCACTGTATGTGTGAGCTATATGGTATTGTGTAATTTTTTTGAAGGAAAATGG
AAATTCCTGAGAACAGTTTGTGTTAAAGAAATATATCAAAATCATTGTGAATAAACTTGATCATCCA
TCTCAATATTGTTTGACATATAAAATAATTATAAGTGTAATAAATTACAATTTAGTGCCACAGTAGTG
AGCATGAAATGAACTATTCAAAAGAGAATATGGCCTGTGCATATTAAAAAATCAAAACAGTGAATGC
AGACTGGAGGAGTAACCTTTGCAAATAAGATGAATATGCTTCATTATAACTCAATATAAAAGGCAAA
TCATCAGAATATTTTAAATGTTGTTGAAAAATGTTTCCCAAGGAAAGTTTATTATTTGCTGCTGTT
TCAAGAAAATTACTTTTACTAAATTTTTTTGTGTGAATTTAAACAGCTAAATAGGGATCAGTAACCTTA
TCTCTATCCTTAATGAACATTTGTTTTATTGGTGGCTGGAATATTCTATTGTATTCTGTGTATATT
TTTAATAAAATTATTTTGGCCTCTTATAAGACAAATCTTATTAATTTGAAACATTTTCATATATACA
CATAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001144903.3
Summary:	May function as sterol isomerase.[UniProtKB/Swiss-Prot Function]
Locus ID:	53346
MW:	30