

Product datasheet for **SC208924**

SAMSN1 (NM_022136) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	SAMSN1
Synonyms:	HACSI; NASH1; SASH2; SH3D6B; SLy2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_022136
Insert Size:	714 bp

Insert Sequence: >SC208924 3'UTR clone of NM_022136
The sequence shown below is from the reference sequence of NM_022136. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AAGATTATTATCACAGAGCCAAGTGACGACACACGCATTCCCAACTATATCTACAGATGCATTCCAT
TTTAACTCTTCTTGAGCTAAAACGTCAAATAGGAGAGGAAGATAAGATAAATATTTGTAATAAACCT
AAAGTTTAAATGTTTTAATCTGAATAATTGTACATAAAATTTGTATCTCTAACATTCCAAATTACTGT
CAATAAAATATATATTTATTATTTAAATGCTATGTGTTAATATTTCACTTGCTTGATTAGAAAGGCA
AAATGTAAGACTTTGGTATGTGTGACATATGCTTTATTTGGCTTTATTTACAAGTACAGTATCTGCAA
AAAACAAAGTAACCTTTTTTCATACCTGCCAGTTTGAATTTATATATGTTATTGAACAAATAGTAATA
GAGGATTCGCTGTTGAAACAAGTTGTCCAAGCAATGTTATATTCATTTTATACTTATTGGGAAAGTGT
GAGTTAATATTGGACACATTTTATCCTGATCCACAGTGGAGTTTTAGTAATTATATTTGTTGATTCT
TCATTTTGTCTTCTGGTATAAAAGTAGAGATAATGTGTAGTCACTCTGATTTAGTGAACCAATTGTA
ATAATTGTGGAAATGTTTTGTCTTTAAGTGTAATATTTTAAATTTGACATACCCTAATGTTAATAAT
AAAAAGAACTATTTGCATATTGCA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTCTATGAAAGG

```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_022136.5</u>
Summary:	SAMSNI is a member of a novel gene family of putative adaptors and scaffold proteins containing SH3 and SAM (sterile alpha motif) domains (Claudio et al., 2001 [PubMed 11536050]).[supplied by OMIM, Mar 2008]
Locus ID:	64092
MW:	27.8