

Product datasheet for **SC208928**

CDK10 (NM_052988) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: CDK10 (NM_052988) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: CDK10
Synonyms: ALSAS; PISSLRE
ACCN: NM_052988
Insert Size: 692 bp
Insert Sequence: >SC208928 3'UTR clone of NM_052988
 The sequence shown below is from the reference sequence of NM_052988. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GAGGGCCAGAGCAAGCGCTGTAACCCTGACGGTGGGCCTGGCACACGCCTGTATCCCACACCAGGTC
TTCCGATCAGTGGTGTCTGTGAAGGGTGCCGCGAGCCAGGCTGACCAGGCGCCGGGATCCAGCTCATC
CCCTTGGCTGGGAACATCCTCCACTGACTTCTCCCACTGTCTGCCCTGAACCCACTGCTGCCCCAGA
AAAAGGCCGGGTGACACCGGGGGCTCCAGCCCGTGCACCCTGGAAGGGCAGGTCTGGCGCTCCATC
CGTGGCTGCAGGGTCTCATGTGGTCTCTCGCTATGTTGAAATGTGCAACCACTGCTTCTGGGAG
GAGTGGTGGTGCAGTCCCCCGCTGTCTTTGAGTTGTGGTGGACGCTGGCCTGGGATGAGAGGGCCCA
GAAGACCTTCGTATCCCCTCTCAGTCGCCCGGGCTGTCCCGTGCATGGGTTGGCTGTGGGGACCCAG
GTGGGCCTGGCAGGACTCCAGATGAGGACAAGAGGACAAGGTATGGGGTGGGAGCCACAATTGAGGAT
ACCCCGAGACTACCAGGAGGCCCTGGGCTGGAGGCTGAGCTGCATCCCTGCTCCCCACATGGAGGACC
CAACAGGAGGCCGTGGCTCTGATGCTGAGCGAAGCTATAGGCTCTGTTGGATAAAAGCTTTTTTAACA
GA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_052988.5</u>
Summary:	The protein encoded by this gene belongs to the CDK subfamily of the Ser/Thr protein kinase family. The CDK subfamily members are highly similar to the gene products of <i>S. cerevisiae</i> <i>cdc28</i> , and <i>S. pombe</i> <i>cdc2</i> , and are known to be essential for cell cycle progression. This kinase has been shown to play a role in cellular proliferation and its function is limited to cell cycle G2-M phase. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2009]
Locus ID:	8558
MW:	24.7