

Product datasheet for **SC208955**

SPINT1 (NM_001032367) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: SPINT1 (NM_001032367) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: SPINT1
Synonyms: HAI; HAI1; MANSC2
ACCN: NM_001032367
Insert Size: 1262 bp
Insert Sequence: >SC208955 3'UTR clone of NM_001032367
 The sequence shown below is from the reference sequence of NM_001032367. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GTCTATAACCACACCACGCGGCCCTCTGAGCCTGGGTCTCACCCGGTCTCACCTGGCCCTGCTTCCTG
CTTGCCAAGGCAGAGCCCTGGGCTGGGAAAACTTTGGAACCAGACTTTGCCTGTTTCCAGGCCAC
TGTGCCTCAGAGACCAGGCTCCAGCCCTCTTGGAGAAGTCTCAGCTAAGCTCAGTCTGAGAAAGC
TCAAAGTTTTGGAAGGAGCAGAAAAACCCTTGGCCAGAAGTACCAGACTAGATGGACCTGCCTGCATAG
GAGTTTGGAGGAAGTTGGAGTTTTGTTTCTCTGTTCAAAGCTGCCTGTCCCTACCCCATGGTGCTAGG
AAGAGGAGTGGGGTGGTGTGACACCCTGGAGGCCCAACCCTGTCTCCCGAGCTCCTTCCATGCTG
TGCGCCAGGGCTGGGAGGAAGGACTTCCCTGTGTAGTTTGTGCTGTAAGAGTTGCTTTTTGTTTATT
TAATGCTGTGGCATGGGTGAAGAGGAGGGGAAGAGGCCTGTTGGCCTCTCTATCCTCTCTTCTCTTC
CCCCAAGATTGAGCTCTCTGCCCTTGATCAGCCCCACCCTGGCCTAGACCAGCAGACAGAGCCAGGAGA
AGCTCAGCTGCATTCCGAGCCCCCACCCTCAAGTTCTCCAACATCACAGCCAGCCCGCCACTGGG
TAATAAAAGTGGTTTGTGGAGTTTCTGGCTCTGTTCTGTGGGCCAGGAGAGGGACTAGGTTGCTGGGG
GTGATGGGACCGTCTGTTCAGCTGGTTCAGGCTTCAGCCACCCTGTTGAGCAGGTTTTCCCATCGTCT
CTCTCGAAACTGGCCTGGCCCTGACTTTGGAGAGCAGTGCCATGTTTGACCACCAGAGGGGAGGAGTG
AGCGGGCCACGGGAGCCACGAGTGTGGGAGGCTGCAGTGTGCTGTGTAGCATTCTTCTGGTGTGA
GCCTCTCCACAAGCCTGTTGAGTGGTTGAGCTCTCCACAATTCTGCTCATTGGAGGTTGGGGTTAGGG
GTGCACAGGCACCTTGGCTCCTACCTCCATGTTAGGGCCAGGTGACTAGCTGTCTCCCTCCCAACAGA
AGCCCTGGCTTCTGCGCTGGAGTGCAGCACCCCAACCAGGCTCTGAGCTCCTTCCAGGAGCCTTTAG
TAGAGACAGCCACTCCAACCCACACCCTGCATTTGTTCAAGTACCTCCCGGTGCCAAAAAAAAAAAA
AAAATCATACTTTGAGCTCA
ACGCGT AAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online >](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001032367.2
Summary:	The protein encoded by this gene is a member of the Kunitz family of serine protease inhibitors. The protein is a potent inhibitor specific for HGF activator and is thought to be involved in the regulation of the proteolytic activation of HGF in injured tissues. Alternative splicing results in multiple variants encoding different isoforms. [provided by RefSeq, Jul 2008]
Locus ID:	6692
MW:	45.2