

## Product datasheet for **SC208958**

### hnRNP R (HNRNPR) (NM\_005826) Human 3' UTR Clone

#### Product data:

|                           |                                                 |
|---------------------------|-------------------------------------------------|
| Product Type:             | 3' UTR Clones                                   |
| Product Name:             | hnRNP R (HNRNPR) (NM_005826) Human 3' UTR Clone |
| Symbol:                   | hnRNP R                                         |
| Synonyms:                 | hnRNP-R; HNRPR                                  |
| Mammalian Cell Selection: | Neomycin                                        |
| Vector:                   | pMirTarget (PS100062)                           |
| ACCN:                     | NM_005826                                       |
| Insert Size:              | 2000 bp                                         |



[View online »](#)

**Insert Sequence:** >SC208958 3'UTR clone of NM\_005826  
The sequence shown below is from the reference sequence of NM\_005826. The complete sequence of this clone may contain minor differences, such as SNPs.  
**Blue**=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
CAGGATACCTTATGGGCAACAGTGAAGAGACAAGTAAGGGCTTGAAAATGATACTGGCAAGATACGAT
TGGCTCTAGATCTACATTCTTCAAAAAAAAAAATTGGCTTAAGTGTTCATCTTTAAGTAGCATTTTGC
TGCCATTTGTATTGGGCTGAAGAAATCACTATTGTGTATATACTCAAGTCTTTTATTTTCTCTTTT
CATAAATGCTCTTGGACATTATTGGGCTTGACAGTTCCTTATTCTGGGATTACAATGCTTTTATCG
TTTCAGGCTTCATTTTAGCTTCAAAACAAGCTGGGCACACTGTAAATCATGATTTTGCAGAACCTTTG
GTTTTGGACAGTTTCATTTTTTTGGATTGGGATAGATTACATAGGAGTATGGAGTATGCTGTAATAA
AAATACAAGCTAGTGCTTTGTCTTAGTAGTTTTAAGAAATTAAGCAAACAAATTTAAGTTTCTTGTA
TTGAAATAACCTATGATTGTATGTTTTGCATTCTAGAAGTAGTTAACTGTGTTTTTAAATTGTTAT
AACTTCACACCTTTTGAATCTGCCCTACAAATTTGTTTGGCTTAAACGTCAAAGCCGTGACAAATT
TGTCTTTGATGTGATTGTATTTCCAATTTCTGTTTCATGTAAGATTTCAATAAACTAAAAATCTAT
TCAAAACATTACCTATTTCAAATTCATTTGTGCTTAAACATTATTTTATTCGAATCCTTGCAAGGA
ATATTGTTTTCAAAGTATAACTGCCTATGTGAGTGATCAAATTCATGCAAAATTTCTGTTTATCCAA
CATGTGTTGTGGATATTAGACTGAGGATTAGTCTAAATTCAGATTTTATTCGAATAGCACTTTTCATT
GATTGGTGTCTCTCATAGTCTCTTCTAATTTTGAATTTTATCTTTTATTAACCAATAGCACTTTTCATT
TCATAGACTTTTGAGCATTCTTGACGGCATTGGATAATTTGCATTAACAAATAAAATGTGGCTGGGA
AAGCAGACTGCAATTTATTTTGACCCACATTCTTTATAAGGGATAATTTGTGACTGAAGACCTTATTTG
CTGGATAAGTGAATGTACATACAGGGAATTTGTTAGTACTGTGGCGTCAATTGCACTGCGGACTTTAG
TTTTAGTGAAATGATCGTGACTACCATTTCACTTACTGGTTAGGAATATAATTTCTATTACTATGT
AAATTTGCTTTGAGAACATTTCATGCCTTGTAAGAGAGGCTTAAAAATATCCAGTAGTAATACCTTTGA
CTTTCTGTTAACACTTCGTTTATTTTATTTATCTGCTATCCCTGATGCCCAATTTTGTGGCCAA
GACATAACATCTATTTCCCTATTTCTTTAGAAATCTGATTTATTCTAATGACTTTTTTCTTGAGAA
TGAGAGGAACAAAGAATCTAAGCTTTACGGTACAGAAAAGTGAATCACTTTCAAATTACTATTAAGTT
TTCCCTCCATCTTATTCTACGTGGTTCTCAAATGTTTTTATAATGATAAAATATCAGAAGCCATTG
GGGAATAATTAATTAATGGATACACCATTTCTGTGTTTTGTTCTGTTTGGTGTCTTTGTATAGGAAT
ATTAGATTTCAAATCATAAGCCACTTAATCCGTGTAATCTTGGCTAGATCATATTGGGTTTCATGGC
TTCTTTTCATTTCAATTTAAGATTACTATTTACTGATGAAAGTCAATATATTTTCATTTCTTTTCAAAG
TACACTCATTACACAGTATAGTCATTTAACAGACCTTATTGATATTTTAGTTTTTGGAACTGATAAA
GGAATTATGGTCTGTTGTATGATAGTTTTGAGTTAACATGGGTATTTCTTTGGTAAGACAGATCCAGTT
GTTTTTCTTATGCCCATGTATGCAAGTAAAGACACCAACATTCTTTAAATGCTGCTTTTAACTTTG
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

**Restriction Sites:** Sgfl-MluI

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

**Components:** The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

**RefSeq:** [NM\\_005826.5](#)

|           |                                                                                                                                                                                                                                                                                                                                                                                 |
|-----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Summary:  | This gene encodes an RNA-binding protein that is a member of the spliceosome C complex, which functions in pre-mRNA processing and transport. The encoded protein also promotes transcription at the c-fos gene. Alternative splicing results in multiple transcript variants. There are pseudogenes for this gene on chromosomes 4, 11, and 10. [provided by RefSeq, Jul 2014] |
| Locus ID: | 10236                                                                                                                                                                                                                                                                                                                                                                           |
| MW:       | 78.2                                                                                                                                                                                                                                                                                                                                                                            |