

Product datasheet for SC208962

C16orf71 (NM_139170) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: C16orf71

Mammalian Cell Selection: Neomycin

Vector: pMirTarget (PS100062)

ACCN: NM_139170

Insert Size: 711 bp

Insert Sequence: >SC208962 3'UTR clone of NM_139170
The sequence shown below is from the reference sequence of NM_139170. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GCCCTGATGCCTCCTCTGGAGCAACTATAGCTGCCTCAGGATGTGCTGCTGTCTGCCCGTGGAAG
AGCAGAGAAATCATCACCACTTGGGCCCCACGGTGCCACGGGCTCAGGCAGCACAGTAGGGCGCCGGG
CCTTTGGGACAGTGTCCAGCTTCCCTGGGGTTCACCCCTGGCTGCCAGGCCACTGAGGATGGGCAGG
GGTCTCTTCTCATCAAGCCTTGTACCAGGCAAAAGACAGGCCCTGCTTGGCCGTGGTCACTGGCCGCCA
AGATCAGGGCTACAGATGTCTGCTCTCTGGACCCACGTGATCTGGCCACTGGGACCCCAACCCGACC
CCACTCCCAATGATGAGGGGCATTTTCATTGCAAGTCAAAGGCAAGACAGGCTTCCATAAAGTCCAGG
AGGTCCCTCCTGTAGGGCACAAGGCCAGGCTGCCTCCAGCCCCAGGCCCTCTCCACCTACAGAGAC
CCTCCCTGCCCTCCACTCCGGGGCCTGTGCCGCCAGAACGGGCTCTGCCCCATACGCTGCCCTC
GCAGCCTGGCGGCCTCCGCTGTGGCTGCCTAGCTGTCAAGAGCAAAGGCTTTTTTTTCTTCAACCCCA
TTTTCTCCATTTCTCCACCTTTTAATGCCAGTAACCTCACTGAGAATGTTTACAGTGATGAAAA
TAAACTCTGTTCCAAGTTCAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_139170.3
Summary:	In ciliated cells, dynein axonemal particle-specific protein required for deployment of ODA to the axoneme. Interacts with outer dynein arm (ODA) subunits.[UniProtKB/Swiss-Prot Function]
Locus ID:	146562
MW:	26