

Product datasheet for **SC208970**

SH3YL1 (NM_015677) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: SH3YL1

Synonyms: RAY

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_015677

Insert Size: 697 bp

Insert Sequence: >SC208970 3'UTR clone of NM_015677
The sequence shown below is from the reference sequence of NM_015677. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
TTTCCAGCCAACACTACGTAACCATGAATAAAGCGTATACTATTTCTTCTTTGAGAATTACAAAAAAT
TATTTCTACACTGACAGGATTTACTAGTTAAGCAATGTTTAAATATAAATTTAAAAAATCTGTCTTA
CAAAATTTCCATTCCGTATGTAAGAGATTTGTTTTCTATATAAAAGAGCTGACTGACATATCTTTA
AATACTTTGTACTAATTTATCACACTTACTGTGTCATAGAATATCATACAGTTTATACGCTCATAGTT
CTCTTGTAACACTTCAAACATCGCTAAGCATTGATCTGGCCATGTATATGGTAGCTGTGTTTAATT
TGAGAACTTGAGGGTAGAGCCACAAATTTCAATTCCTACATTTCCATTTGCAAAGTGACTAGAGAAAA
AGAAATCAGCTTAAATGAGGTATTAAGTAATGTTTAGAGTCGTAGGTATTAAGTAGAATATAAATCCTT
AGAAATTGTCTTTATACCTTCAAAATTTACTATGCATTTATCATAGAAATGTGATTACAAAGAAAGTC
TGACTACCATGTCTTTAAACATATGGCATCTCTCAACTTTCTTCTTATGGGGCTACATTTGTTTCATT
TCCAGCAGTAGCATAACTTACGGTGACATGGTAGACTGTCTCTAAATAAAATTTTAAATGTTTACT
AAGTTAA
ACGCGTAAGCGGCCCGCGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_015677.4</u>
Locus ID:	26751
MW:	27.4