

Product datasheet for **SC208981**

FAM149A (NM_001006655) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	FAM149A
Synonyms:	MST1I9; MSTP1I9
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_001006655
Insert Size:	723 bp
Insert Sequence:	>SC208981 3'UTR clone of NM_001006655 The sequence shown below is from the reference sequence of NM_001006655. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC TTCAAGAGGAGATTCCAAGTGACATCTGAACCACCTCACCAGAACCATTGGAGCACCTGTGGCCTCG GCACACTGCTTATCACTTGAAAAATGGAGGAACAAGTGTTATTTATCTGTGTCTGACAGTGTGAG ATGTTAGACCGAGAGAAAAGCAAACAATAAATCACTTAATCTTGAACTTTGCACTGTAAGAGAGT GAAAGTCAAACCCACCAAGCTGCCTTGCTGAAAGCATCTCCAAGTTCTGTAGGCTTTGTTCAAGCC CCTGATGTGTGTTGCAAGTCCGTCATTTTATCAATGATACGCAACATAATTAGCATGTGTATGTA ACACACTAATTACATGAAGAATATATTATCTAGGAATATGTCATTTGTAACCTTTGACTCTGACGTAT TTTCTAGTCTCACCCTATTAGTAATTCTGTCAAAATCCTAGACGCTAGTAGATATTAGTATTTGAGACA ATATGTCAGAGTCTGAATAAGACTCAGCTTGTTCTAGTCATGCTGAGTTTCTCATATATTTTTTAAAG TGCTTCCGCCATGAACCTCTTTGAGAAATGTCAGTGGCAAAATCTCCCAAGCTAATTTTCTTTTCTAT GTTCCATCTTCAAATTCAGACAGTTCACGTGAAGTACCCACATTTATTGATCTTAGGGTTGAATATA ACGCAATTAAGCTAAAAAAAAAAAAAAAAAAAA ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001006655.3</u>
Locus ID:	25854
MW:	28.1