

Product datasheet for **SC208985**

ZNF626 (NM_145297) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF626 (NM_145297) Human 3' UTR Clone
Symbol:	ZNF626
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_145297
Insert Size:	727 bp
Insert Sequence:	>SC208985 3'UTR clone of NM_145297 The sequence shown below is from the reference sequence of NM_145297. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
GAGAGCCAAAGTCCTCTTCATGACATAAAGAGACTGCACAGTGTGGCTGCTTTCCATAATTTTGGG
ACGCACAAATGTCTGCATGATTTTGAGAACTAAATTTACTCTCAAAGTTCTCTTGCATCAGATCTGA
AATGTCTGAGAGTAATAGTTTCTGTTGAATTTTTTTTGTTCATTTTCTGCACAGTCCATTCTGTTT
TATTACTATCTAGGCTTGAAATATATAGTTTGAAATTATGACATCCTTCCTCTTTGTTATTTTCCATCAT
GATTGCTTTGGCTATTCAAAGTTTATTTTAGTTTCATGTAATTTTGAATTGTATTTTCCATTATTGT
GAAAATAGTACCACTGCAATTTTAATAGGAAGTTTATTGAATCTATAGATTACTTTGGATAATATGGCA
CTTCAATAATATTCATGTTTTCAATTCATAGACAAAATATTTTAAATTTATTTGTATCTTTTCTAACT
TTTCCTTTTTTTTATTGTAAAGATTTACCTCCTTGGTTAATATTTTCTCAGAAATTTATTATTTAAGGT
ATAGTCAATAAAATTTTCTCTCTATTTTGTGATAGTTTAAAGTGTATGAAACCATAGATATACTTG
TATGTTAATTTTATTTTGTCTAATTTACTGAGTGTATTTATTAGTTTAGAGAGGTTTTAATGTACTGT
TTATGGTTTTTTTAAATATAAGATTACTTATTTTTTAA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_145297.4
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	199777
MW:	29.2