

Product datasheet for **SC208992**

C17orf65 (NM_178542) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: C17orf65

Synonyms: DKFZp762C2414; MGC119234

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_178542

Insert Size: 724 bp

Insert Sequence: >SC208992 3'UTR clone of NM_178542
The sequence shown below is from the reference sequence of NM_178542. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AAGGCTTGGACCCCTTCAGCTTGCTGTGACACTCCCTTGCCTTTCCACCCCTCTCAGCTTGAGGTCCT
ACCTGGTGTCTCCATGAGAATCATACAAGCCCGCTACAGACAACAGAATTGGAAGGTCCTGAGAGCT
ACCAGAAGCATCTTCAGTTTTATATGAAAACTTAAGCCAGGAGTGATTTGTGCAGCCAAGGACACA
CAGCAAGTTAGTGGCTGCACCAGAAATAGTACTCACTGCCTCAGACAGTCCCTCCACAGTCTACTCTG
TCTCCACCCCAAGGCTCTGGTTCTGAATCATTAGTTGTGAGTTGTCTCTTTTTTTTTTTTTTTTTTTGG
AGACAGGGTCTCTCTGTGGCCAGGCTGAAGTGCAGTGGTGTGATTACAGCTCATGTAGCCTGGACC
TCCCAAGCTCAAGCAATCCTCTGCCTCTGCCTCCCAAGTAGCTGGGACCACAGGTGGGCAGCACCACA
CCTGGCTAATTTTTTAAATTTTTTTTATAGAGACAGGCTCTTACATTTGCCCCAGCTGGTCTCAAAC
CCCTGGGCTCAAGCAAGCTCTGCCTTGGCCTCCCAAGTGTGGGATTACAGGTGTGAGCCACCACA
CCCGGCCTTATCTCTAACATTTTGATAGCACAAAATTTGAAAACCCACTTCTGTGTACGTTATGAAG
CATAGAAAAATTAACAGTATGCACCAAAAAAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTCTATGAAAG
```

Restriction Sites: SgfI-MluI



OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_178542.3</u>
Locus ID:	339201
MW:	27.1