

Product datasheet for **SC209039**

CABYR (NM_012189) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: CABYR (NM_012189) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: CABYR
Synonyms: CABYRa; CABYRc; CABYRc/d; CABYRe; CBP86; CT88; FSP-2; FSP2
ACCN: NM_012189
Insert Size: 729 bp
Insert Sequence: >SC209039 3'UTR clone of NM_012189
The sequence shown below is from the reference sequence of NM_012189. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GAACCAGAAGGGGAATCAACAGCTGAATAGGTTTGATGAAGCCAGCAATGGCAACAAGTGAACGAGGA
CAACCACCACCATGTTCTAACATGTGGACCTTTATTGTCTAACTGATAAGAATCAACAAGTCCACCA
TCACCGCCACCTGCACCTGGGCCTTTTCCCAAGCAACCCTCTATTTACCTAATCCTAAGGATCCACAG
TTTCAGCAGCATCCACCAAAAGTCACTTTTCCAATTATGTGATGGGCGACACCAAGAAGACCAGTGCC
CCACCTTTTATCTTAGTAGGCTCAAATGTTTCAGGAAGCACAGGGATGGAAACCTCTTCTGGACATGCT
GTCGTTTCACAGTCAGATGCTTGAGATATGTTGCAATGCAAGTGCCATTGCTGTTCTGCAGATGAG
AAATACCAGAAACATACCCTAAGTCCCCAGAATGCTAATCCTCCAAGTGGACAAGATGTCCCAGGCCA
AAAAGCCCTGTTTCTTTCTGTTGCTTTCCAGTAGAAGATGTAGCTAAAAAAGTTCAGGATCTGGT
GACAAATGTGCTCCCTTTGGAAGTTACGGTATTGCTGGGGAGGTAACCGTGACTACTGCTCACAAACGT
CGCAAAGCAGAACTGAAAAGTATCCAGAAATGACGCTGTCTGGGTCAACATTTTCAGGGAGGAGTCTG
CCACCAGTGAATGTATCAATAAACTTCATGCAAGCATA
ACGCGTAAGCGGCCGCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites:

Sgfl-Mlul

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_012189.4</u>
Summary:	To reach fertilization competence, spermatozoa undergo a series of morphological and molecular maturational processes, termed capacitation, involving protein tyrosine phosphorylation and increased intracellular calcium. The protein encoded by this gene localizes to the principal piece of the sperm flagellum in association with the fibrous sheath and exhibits calcium-binding when phosphorylated during capacitation. A pseudogene on chromosome 3 has been identified for this gene. Alternatively spliced transcript variants encoding distinct protein isoforms have been found for this gene. [provided by RefSeq, Jul 2013]
Locus ID:	26256
MW:	26.2