

Product datasheet for **SC209045**

C11orf52 (NM_080659) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: C11orf52

Mammalian Cell Selection: Neomycin

Vector: pMirTarget (PS100062)

ACCN: NM_080659

Insert Size: 702 bp

Insert Sequence: >SC209045 3'UTR clone of NM_080659
The sequence shown below is from the reference sequence of NM_080659. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TATGACAGCAAGAACGGGACCTGGTGTGAGCGCTTGGGAGGAAGGCCAGTCCATCGTTAACCCTAC
ACCTGTGGGGGAGAACCTACTGCTTTGGGGAATTGGGTGGCAACCCAGGGATGTTGTCCACGTTTTAAG
CTTAAAGAACTCCAAAGCCCTGGAATTGTGGGCGTCCCATTTTCTCCCTGGCCTCTTACTTGCCACT
GTTAGGTTGGAGTTAATAGTTGGTTAGCTATGGGTGGATAGCTAACTTAGCTATCCACATGAGGTTA
GGTGGAGTGTGCAGGGAGGTAGGTCTTCGACCCACCATTTGCTCCTGCTCTTGCAAAGTGCAGAAA
GAGCAGGTAAAAAGTTAGCCCATAGCACGTGGGCTGATCTTGTTGGACTTTAATTAATGGTATCCTTT
TTCACACACCTTAACTCCAAAGCTGAAGGGGAGGTGCTCCGAAAAGCATTCTGAGCTTTTACATTCA
GACACTAGGGCTGGCTTGGCCATTTCTCAGTGGCTGAGTCAGGTTCCAGACCTTAACAGACAAATGG
GGGAAGGGGCTGAGGCGTGGAGGAAACGCTTTCAGGACCTCAAGTCAGTGATAACCGACAGCACAGTG
AGGCCTGGAAACTCCCTCGGAGGACTGTTGGCTCTCGCAGGTGCGGGGCCGCTAGGATCAATAAACG
GCCTGAGAAACA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_080659.3</u>
Locus ID:	91894
MW:	25.9